

ApexPrep DNA Plasmid Miniprep Kit (50 Tests)

Please add 200 μ L RNase A to Buffer A1 before use, then store Buffer A1 at 2-8°C.

Please add 48 mL ethanol to Buffer AP before use.

1. Inoculate 1-5 mL of LB medium with selected antibody added with one single colony, culture at 37°C overnight with vigorous shaking.
2. Collect cells by centrifuging at 3000g for 5 minutes. Discard supernatant and keep the pellet.
3. Resuspend pellet in 250 μ L Buffer A1. Transfer the mixture to a 1.5 mL microcentrifuge tube.
4. Add 250 μ L Buffer A2 to the mixture. Invert 4-6 times to mix well.
5. Add 350 μ L Buffer A3 to the mixture. Invert 4-6 times to mix well.
6. Centrifuge at 10,000g for 10 minutes. Insert a spin column to a collection tube. Transfer supernatant to the spin column and discard pellet.
7. Centrifuge spin column with collection tube for 30 seconds. **Transfer flow through to the spin column and centrifuge again.**
8. Add 750 μ L Buffer AP to the column and spin for 30 seconds. Discard flow through.
9. Centrifuge for another minute to remove all residue buffer.
10. Transfer the spin column to a clean 1.5 mL microcentrifuge tube. Add 30 μ L Buffer AE to the center of the membrane. Let it stand for 1 minute.
11. Centrifuge for 30 seconds and collect flow through. This is your eluted plasmid DNA.