

Product Information

NADH Fluorometric Assay Kit

I. Kit Contents:

Components	K2037-100 100 assays	Cap Color	Part Number
NADH Extraction Buffer	50 ml	NM	K2037-C-1
NADH Cycling Buffer	15 ml	NM	K2037-C-2
PicoProbe™ (DMSO)	0.4 ml	Blue	K2037-C-3
NADH Cycling Enzyme Mix (Lyophilized)	1 vial	Green	K2037-C-4
NADH Standard (Lyophilized)	1 vial	Yellow	K2037-C-5

II. Introduction:

Nicotinamide adenine dinucleotide (NAD) is a coenzyme present in all living cells and exists in two forms: NAD⁺ and NADH. NAD plays important roles in energy transforming and redox state of cells or tissues. NADH, the reduced form of NAD, plays an important role in cell regulation and repair processes. Measurement of low level NADH in samples or in enzymatic reactions is of increasing interest.

The NADH Fluorometric Assay Kit provides a highly sensitive and convenient way for detection of low level NADH in samples or in enzymatic reactions based on fluorescence method. The NADH Recycling Enzyme Mix specifically recognizes NADH in the enzyme recycling reaction, which are not required to purify from samples. The assay is fast, simple and convenient, and can measure less than 8 nM NADH in a variety of samples.

III. Application:

Measurement of NADH in various tissues/cells.

Measurement of low levels of NADH in enzymatic reactions.

Analysis of metabolism in various cells.

IV. Sample Type:

Animal tissues: Liver, muscle, heart etc.

Cell culture: Adherent or suspension cells.

Enzymatic reactions.

V. User Supplied Reagents and Equipment:

96-well white plate.

Multi-well spectrophotometer (Fluorescence reader).

VI. Storage and Handling:

Store kit at -20°C, protected from light. Warm NADH Extraction Buffer & NADH Cycling Buffer to room temperature before use. Briefly centrifuge small vials prior to opening.

VII. Reagent Preparation and Storage Conditions:

PicoProbe™: Ready to use as supplied. Warm to room temperature before use. Store at -20°C.

NADH Cycling Enzyme Mix: Reconstitute with 220 μ l NADH Cycling Buffer. Pipette up and down to dissolve completely. Aliquot and store at -70°C . Avoid repeated freeze/thaw. Keep on ice while in use.

NADH Standard: Reconstitute with 200 μ l DMSO to generate 1 mM (1 nmol/ μ l) NADH Standard solution. Keep on ice while in use. Store at -20°C . Use within two months.

VIII. NADH Assay Protocol:

1. Sample Preparation: Liquid samples can be measured directly. Tissue (~ 10 mg) or cells ($\sim 1 \times 10^6$) should be rapidly homogenized with 200 μ l ice cold NADH Extraction Buffer for 10 minutes on ice. Centrifuge at 12000 rpm for 5 min. Collect the supernatant. Add 1-50 μ l samples into an eppendorf tube/white plate and bring the volume to 80 μ l with NADH Extraction Buffer.

Notes:

- For unknown samples, we suggest testing several doses of your samples to ensure the readings are within the standard curve range.
 - Cell or tissue lysates may contain enzymes that consume NADH rapidly. We suggest removing these enzymes by filtering the samples through 10 Kd molecular weight cut off filters before performing the assay.
- NAD Decomposition: To detect NADH, the NAD needs to be decomposed before the reaction. Put samples at 60°C for 30 minutes to completely decompose the NAD. Cool samples on ice. Centrifuge briefly and transfer 50 μ l of samples into a 96 well white plate.
 - Standard Curve Preparation: Dilute NADH to 10 pmol/ μ l by adding 10 μ l of 1 mM NADH to 990 μ l dH_2O , mix well. Dilute further to 0.1 pmol/ μ l by adding 10 μ l of 10 pmol/ μ l NADH into 990 μ l dH_2O , mix well. Add 0, 4, 8, 12, 16 & 20 μ l of diluted 0.1 pmol/ μ l NADH Standard into a 96 well plate to generate 0, 0.4, 0.8, 1.2, 1.6, and 2 pmol/well NADH Standards. Adjust the volume to 50 μ l/well with NADH Extraction Buffer.
 - Reaction Mix: Mix enough reagents for the number of assays (samples and standards) to be performed. For each well, prepare 100 μ l Reaction Mix containing:

	Reaction Mix	Background Control Mix
NADH Cycling Buffer	96 μ l	98 μ l
NADH Cycling Enzyme Mix	2 μ l	--
PicoProbe™	2 μ l	2 μ l

Add 100 μ l of the Reaction Mix to each well containing the Standard & test samples, mix well.

If your sample has fluorescence background, prepare a parallel sample well as the background control.

5. Measurement: Incubate the reaction for 30 min. at room temperature. Measure fluorescence (Ex/Em = 535/587 nm).

6. Calculation: Subtract the 0 NADH Standard reading from all readings. Plot the NADH Standard curve. If the sample background control reading is significant, subtract the background control reading from the sample. Apply the corrected sample reading to the Standard Curve to get B pmol of NADH in the sample well.

Sample NADH concentration (C) = $B/V \times \text{Dilution Factor} = \text{pmol}/\mu\text{l} = \text{nmol}/\text{ml}$

Where: B is the amount of NADH in the sample well (pmol).

V is the sample volume used in the reaction well (μ l).

NADH in Samples can also be expressed in pmol/mg of sample.

NADH molecular weight: 663.43 g/mole.




Figure: NADH Standard Curve (a). Decomposition of NAD, but not NADH (b) & measurement of NADH in Rat kidney (1.74 μ g), muscle, (2.51 μ g), heart (1.59 μ g) and liver lysates (29.5 μ g) (c).

For research use only! Not to be used in humans.

Our promise

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For more details, please visit <http://www.apexbt.com/> or contact our technical team.

Tel: +1-(832)696-8203

Fax: +1-832-641-3177

Email: sales@apexbt.com