

## Product Information

### BCA Protein Quantitation Kit

#### I. Kit Contents:

Components	K2185-1000 1000 assays	Cap Color	Part Number
BCA Reagent	100 ml	NM	K2185-C-1
Copper Reagent	2 ml	Blue	K2185-C-2
Bovine Serum Albumin (10 mg/ml)	1 ml	Yellow	B7773

#### II. Introduction:

The BCA Protein Assay Kit offers a fast, convenient and detergent tolerant method for measuring the concentrations of proteins in solution. A copper ( $\text{Cu}^{2+}$ ) salt is reduced to cuprous state by protein in the BCA Protein Assay Kit. The produced  $\text{Cu}^{2+}$  ion forms a strongly colored complex with the bicinchoninic acid reagent showing an intensive absorbance band centered at 562 nm. The intensity of the blue complex is proportional to the amount of protein in the sample. The concentration range of measurement for this kit is 0.5 - 30  $\mu\text{g}$  protein (0.01 - 0.6 mg/ml).

#### III. Storage Conditions:

The BCA and Copper Reagents are stable at room temperature. The Bovine Serum Albumin should be aliquoted after the first thaw and stored at  $-20^{\circ}\text{C}$ . All reagents are stable for up to 12 months under proper storage conditions.

#### IV. Reagent Preparation:

Prepare Working Solution by adding 1 part of Copper Reagent to 50 parts of BCA Reagent. The total volume made will depend upon the number of samples and standards to be quantitated. Each sample and standard will require 100  $\mu\text{l}$  of working reagent. Once made, the Working Solution is stable for several days.

#### V. Protein Assay Procedure:

The BCA protocol is very flexible. Both the incubation time and temperature can be varied over a rather wide range. Lower protein samples can be more easily quantified using higher temperatures and longer incubation times.

1. Preparation of Bovine Serum Albumin Curve: Label (8) tubes 1 - 8. Dilute the Bovine Serum Albumin to 1 mg/ml Stock Solution (i.e., 40  $\mu\text{l}$  + 360  $\mu\text{l}$  buffer). Ideally, use the same buffer contained in your samples. Transfer 200  $\mu\text{l}$  from tube 8 to tube 7. Then prepare serial dilutions as below:

Tube Code	BSA solution ( $\mu\text{l}$ )	Buffer, $\mu\text{l}$	50 $\mu\text{l}$ =
8	Stock Solution (256)	144	32 $\mu\text{g}$
7	Tube 8 (200)	200	16 $\mu\text{g}$
6	Tube 7 (200)	200	8 $\mu\text{g}$
5	Tube 6 (200)	200	4 $\mu\text{g}$
4	Tube 5 (200)	200	2 $\mu\text{g}$
3	Tube 4 (200)	200	1 $\mu\text{g}$
2	Tube 3 (200)	200	0.5 $\mu\text{g}$
1	---	200	0 $\mu\text{g}$

2. Dilute samples to fall within 0.01 - 0.6 mg/ml range.

3. Pipette 50  $\mu$ l Standards or samples into duplicate wells in a clear bottom 96 well plate.
4. Add 100  $\mu$ l of Working Solution into each well that contains the standard or samples.
5. Shake gently to mix. Incubate for between 30 - 90 min at 37°C - 60°C. Cool to room temperature.
6. Measure OD at 562 nm. The signal is stable for at least 1 hour.


## VI. Calculation:

Subtract the blank OD (0 standard) from all standard and sample OD values. Plot the corrected OD against standard protein concentrations. Use the standard curve to determine the sample protein concentration. Alternatively, the equation for the best line fitting the standards can be used to determine the protein concentration of your samples.

## VII. General Considerations:

For unknown samples, several dilutions of a sample should be tested to ensure the OD reading is within the standard curve range.

When assaying protein in solutions containing detergent, best results are obtained by adding the same amount of detergent to the wells containing the protein standard.



**For research use only! Not to be used in humans.**

## Our promise

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For more details, please visit <http://www.apexbt.com/> or contact our technical team.

Tel: +1-(832)696-8203  
 Fax: +1-832-641-3177  
 Email: [sales@apexbt.com](mailto:sales@apexbt.com)