

CFDA SE Cell Tracer Kit

1. Introduction

The CFDA SE Cell Tracer Kit provides a versatile and well-retained cell tracing method. The kit is simple and convenient for users which containing CFDA SE dye and DMSO.

CFDA SE can passively diffuse into cells. It is colorless and nonfluorescent until entering a live cell, the acetate groups of CFDA-SE are cleaved by intracellular esterase to create amine-reactive and highly fluorescent carboxyfluorescein succinimidyl ester (CFSE). The succinimidyl ester group would react with intracellular amines when entering into cells, forming fluorescent conjugates that are well-retained and can be fixed with aldehyde fixatives. Redundant unconjugated reagent and by-products passively diffuse to the extracellular medium, where they will be washed away.

The dye-protein adducts that formed in labeled cells will be retained by the cells throughout development, meiosis, and in vivo tracing. The label can be inherited by daughter cells after cell division or cell fusion. It isn't transferred to the adjacent cells in one population. Cell division will result in sequential halving of CFDA SE fluorescence appeared in a cellular fluorescence histogram in which the peaks represent successive generations. In previous literature, viable hepatocytes labeled were easily located by fluorescence microscopy even 20 days after intrahepatic transplantation.

CFDA SE Cell Tracer Kit contains dry dye and DMSO. Users should be careful to extend the shelf life of the dye by not dissolved in solution. A stock solution is prepared by dissolving the dye in anhydrous DMSO until prior to use.

2. Experimental Protocol

2.1. Guideline

The following protocols describe introducing the CFDA SE into the cultured cells and imaging the stained cells by fluorescence microscopy. Our recommended initial conditions might require modifications because of differences in cell types, culturing conditions, desired effect, and so on. The optimal concentration of dye for staining is variable depending upon the application; we recommend preliminary experiment in a tenfold range concentration gradient. In general, long-term staining (more than about three days) or the use of rapidly dividing cells requires 5–10 μ M concentration. Lower concentration (0.5–5 μ M) is needed for shorter experiments, such as viability assays.

To maintain normal cellular physiology and reduce potential artifacts from overdose, the concentration of dye should be kept as low possible.

Note: The CFDA SE dye would react with amine groups, therefore, don't use with amine containing buffers or lysine-coated slides.

2.2. Preparing the Reagent

Before opening the vial, allow it warm to room temperature. Prepare an approximately 10 mM CFDA SE stock solution immediately prior to use by dissolving 1 mg CFDA SE dye in 179 μ L high-quality DMSO provided. Dilute the stock solution in phosphate-buffered saline (PBS) or other suitable buffer to the desired working concentration (0.5–25 μ M). Be careful that solutions of the reagent should be used immediately.

2.3. Labeling Adherent Cells

- Grow cells on coverslips inside a petri dish filled with the appropriate culture medium.
- When the cells have reached the desired density, remove the medium from the dish and add preheated (37°C) PBS containing the dye(prepared in step 2.2)..
- Incubate the cells for 15 min at 37°C.
- Replace the loading solution with fresh, preheated medium and incubate the cultures for another 30 min at 37°C.

Tel: 021-55669583; Fax: 021-55669583

Web: <http://www.apexbt.com/>; Email: sales@apexbio.cn

During this time, CFDA SE will undergo acetate hydrolysis. If the cells are to be fixed and permeabilized, continue to step 2.5.

2.4. Labeling Cells in Suspension

- a. Centrifuge to obtain a cell pellet and aspirate the supernatant.
- b. Resuspend the cells gently in preheated (37°C) PBS containing the dye (prepared in step 2.2).
- c. Incubate the cells for 15 minutes at 37°C.
- d. Re-pellet the cells by centrifugation and resuspend in fresh preheated medium.
- e. Incubate the cells for another 30 minutes at 37°C to ensure complete modification of the dye and then wash the cells again. If the cells are to be fixed and permeabilized, continue step 2.5.

2.5. Fixing and Permeabilizing

- a. Before fixation, the cells must be washed with PBS or other suitable buffer.
- b. Standard fixation protocols using aldehyde-containing fixatives should effectively crosslink the amines of the protein–probe conjugate. Typically, we fix the cells for 15 minutes at room temperature using 3.7% formaldehyde.

2.6. Flow cytometry

For researchers who wish to analyze labeled cells and/or study cell division via flow cytometry, we recommend the excellent protocol described in reference [7], [8].

3. References

- [1]. Bronner-Fraser M. Alterations in neural crest migration by a monoclonal antibody that affects cell adhesion. *J Cell Biol.* 1985 Aug;101(2):610-7.
- [2]. Nose A, Takeichi M. A novel cadherin cell adhesion molecule: its expression patterns associated with implantation and organogenesis of mouse embryos. *J Cell Biol.* 1986 Dec;103(6 Pt 2):2649-58.
- [3]. Lyons AB1, Parish CR. Determination of lymphocyte division by flow cytometry. *J Immunol Methods.* 1994 May 2;171(1):131-7.
- [4]. Hodgkin PD1, Lee JH, Lyons AB. B cell differentiation and isotype switching is related to division cycle number. *J Exp Med.* 1996 Jul 1;184(1):277-81.
- [5]. Weston SA1, Parish CR. New fluorescent dyes for lymphocyte migration studies. Analysis by flow cytometry and fluorescence microscopy. *J Immunol Methods.* 1990 Oct 4;133(1):87-97.
- [6]. Karrer FM1, Reitz BL, Hao L, Lafferty KJ. Fluorescein labeling of murine hepatocytes for identification after intrahepatic transplantation. *Transplant Proc.* 1992 Dec;24(6):2820-1.
- [7]. Hou WS1, Van Parijs L. A Bcl-2-dependent molecular timer regulates the lifespan and immunogenicity of dendritic cells. *Nat Immunol.* 2004 Jun;5(6):583-9. Epub 2004 May 9.
- [8]. Dumitriu IE1, Mohr W, Kolowos W, et.al. 5,6-carboxyfluorescein diacetate succinimidyl ester-labeled apoptotic and necrotic as well as detergent-treated cells can be traced in composite cell samples. *Anal Biochem.* 2001 Dec 15;299(2):247-52.