

Direct Mouse Genotyping Kit

Product description

The Direct Mouse Genotyping Kit is specifically designed for species genotype identification. This product contains lysis buffer and balance buffer, which can quickly digest tissues, release complete genomic DNA, and the DNA can be directly used as a PCR template without being extracted from the mixed solution. The 2×Taq PCR Master Mix (With dye) can ensure accurate PCR amplification results. This Master Mix operation only requires the addition of primers, templates and ddH₂O, which greatly simplifies the experimental steps, reduces personal errors and improves the repeatability of the results. This product contains dye, which means it can be directly subjected to electrophoresis after amplification without the need to add sample loading buffer. If agarose Gel electrophoresis of PCR products is required, our product Safe DNA Gel Stain (Item No. A8743) can be added to the gel.

Composition and storage conditions

Components	Size	200 rxns	500 rxns	Storage
2× Taq PCR Master Mix (With dye)		2×1 mL	5×1 mL	-20°C
Balance buffer		20 mL	50 mL	4°C
Lysis buffer		20 mL	50 mL	4°C
Proteinase K		200 μL	500 μL	-20°C
Shipping: Dry Ice		Shelf life: 2 years		

Experimental operation

1. Tissue digestion (taking mouse genotype identification as an example)

- a) Place the mouse tails, toes or ears (~2 mm) in 75 μL Lysis buffer and 0.75 μL Proteinase K for digestion, and incubate at 56°C for 15 min. Then incubate the mixed solution at 95°C for 10-60 min (undissolved tissues do not interfere with PCR).

- b)** After heating, cool the samples to 4°C and add 75 µL of Balance buffer to each sample.
- c)** For every 20 µL of PCR system, use 1 µL of the above-mentioned preparation solution as the template for the next PCR (the sample system can also be adjusted as needed).

***Note:** During the digestion process, shaking the tube several times will contribute to the release of genomic DNA. For most mouse tissue samples, incubation with protease K at 56°C for 15 minutes is sufficient to extract genomic DNA. The tissue may still seem intact, but cracking has occurred. The obtained genomic DNA can be applied to PCR amplification. If not used immediately, the sample can be centrifuged to collect the supernatant and stored at -20°C.

2. Reaction Setup

- a)** Thaw PCR mix, mix thoroughly and then centrifuge. Place on ice for later use.
- b)** Configure the reaction system according to the following table:

Reagent	20 µL Reaction	50 µL Reaction	Final Concentration
2× Taq PCR Master Mix (With dye)	10 µL	25 µL	1×
Forward primer (10 µM)	0.8 µL	2 µL	0.4 µM
Reverse primer (10 µM)	0.8 µL	2 µL	0.4 µM
Template DNA	1 µL	2 µL	
ddH ₂ O	to 20 µL	to 50 µL	

***Note:** Generally, the final primer concentration is recommended to be 0.2-0.4 µmol/L, but it can also be adjusted within the range of 0.1-1.0 µmol/L as needed.

- c)** Invert and mix well, or gently pipette or gently vortex to thoroughly mix the reaction solution, and briefly centrifuge to collect the liquid to the bottom of the tube.

3. PCR Cycling Conditions

- a)** Place the reaction tube on the PCR instrument and start the thermal cycle. Set the reaction program according to the following table:

Procedure	Temperature	Time	Cycles
Initialization	94°C	5 min	1
Denaturation	94°C	30 s	35
Annealing	50-60°C	30 s	
Extension	72°C	30 s/kb	

Final elongation	72°C	10 min	1
------------------	------	--------	---

b) Perform agarose Gel electrophoresis on the PCR products and obtain the results (Our product Safe DNA Gel Stain, item number A8743, is available for your selection).

Notes

Primer Design Guidelines: The primers for qPCR experimental design need to have good amplification effects and few non-specific products. The following design guidelines can be followed:

- 1) The last base at the 3' end of the primer is preferably G or C. The overall distribution of primers A, G, C and T should be as uniform as possible, and areas with high GC or AT content should be avoided.
- 2) The last 8 bases at the 3' end of the primer should avoid consecutive mismatches and the formation of hairpin structures.
- 3) It is preferable that the difference between the T_m values of the forward Primer and the reverse primer does not exceed 1°C, and the T_m value should be adjusted to 55-65° C (It is recommended to calculate the T_m value of the primer using Primer Premier 5).
- 4) Additional sequences of primers, that is, sequences that are not paired with the template, should not be involved in the calculation of primer T_m values. The GC content of the primers is controlled between 40% and 60%.
- 5) Avoid having complementary sequences of more than five bases within the primers or between the two primers, and avoid having complementary sequences of more than three bases at the 3' ends of the two primers.

After the primer design is completed, please use the NCBI BLAST function to search for primer specificity to avoid non-specific amplification.



APEX BIO Technology
www.apexbt.com
 7505 Fannin street, Suite 410, Houston, TX 77054.
 Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com

