

SDS Lysis Buffer

Introduction

SDS Lysis Buffer is a relatively strong cell and tissue lysis buffer. The main components of this product include 50mM Tris (pH 8.1), 1% SDS, as well as various inhibitors such as sodium pyrophosphate, β -glycerophosphate, sodium orthovanadate, sodium fluoride, EDTA, and leupeptin. These effectively prevent protein degradation and maintain native protein-protein interactions. Protein samples obtained using SDS Lysis Buffer can be used for routine PAGE, Western Blotting (WB), and ChIP. This product can be used for animal and plant cell or tissue samples, as well as for fungal or bacterial samples.

Components and Storage

Size	100 mL	Storage
Components		
SDS Lysis Buffer	100 mL	-20°C
Shipping: Blue ice	Shelf life: 12 months	

Protocol

1. Reagents to be prepared by the user

- Protease inhibitors: PMSF (A2587) or Protease Inhibitor Cocktail (K1007)
- Phosphatase inhibitors: Phosphatase Inhibitor Cocktail (K1015)
- Deacetylase inhibitors: Deacetylase Inhibitor Cocktail (K1017)

2. Preparation of lysis buffer

After thawing SDS Lysis Buffer, mix it evenly. Take an appropriate amount of the buffer and add the protease inhibitor PMSF to make the final concentration of PMSF reach 1 mM. Appropriate Cocktails of the above-mentioned protease and phosphatase inhibitors can be added according to the experimental requirements.

***Note:** 1. When more comprehensive protection is required, you can choose to replace it with protease inhibitor Cocktail (K1007);
2. If the protected protein is in phosphorylated or acetylated form, the phosphatase inhibitor Cocktail (K1015) or the deacetylase inhibitor Cocktail (K1017) can be used, respectively.

3. Cell or tissue lysis

3.1 Adherent cells

- a. Discard culture medium. Wash once with PBS, saline, or serum-free medium (washing may be skipped if serum proteins do not interfere).
- b. Add lysis buffer at 150–250 μL per well of a 6-well plate. Pipette gently to ensure full contact.
- c. Most animal cells lyse within 1–2 seconds. Plant cells are recommended to lyse on ice for 2–10 min. For ChIP, after initial lysis, continue lysis for an additional 10 minutes in an ice bath.
- d. After complete lysis, centrifuge at 10,000–14,000 $\times g$ for 3–5 min. Collect the supernatant for downstream PAGE, WB, ChIP, etc.

3.2 Suspension cells

- a. Collect cells by centrifugation. Disperse cells by gentle vortexing or flicking the tube bottom.
- b. Add lysis buffer at 150–250 μL per well-equivalent of a 6-well plate. Flick the tube to aid complete lysis. No obvious cell pellet should remain after lysis. For large cell numbers, split into 0.5–1 $\times 10^6$ cells per tube before lysis. For ChIP, after initial lysis, continue lysis for an additional 10 minutes in an ice bath.
- c. After complete lysis, centrifuge at 10,000–14,000 $\times g$ for 3–5 min. Collect the supernatant for downstream PAGE, WB, ChIP, etc.


3.3 Bacteria or yeast

- a. Take 1 mL culture, centrifuge, and discard supernatant (may wash once with PBS to remove residual liquid). Disperse by gentle vortexing or flicking.
- b. Add 100–200 μL lysis buffer, mix gently, and lyse on ice for 2–10 min. For improved lysis: treat bacteria with lysozyme and yeast with lyticase before using this buffer.
- c. After complete lysis, centrifuge at 10,000–14,000 $\times g$ for 3–5 min. Collect the supernatant for downstream PAGE, WB, ChIP, etc.

3.4 Tissue Samples

- a. Cut tissue into small pieces.
- b. Add lysis buffer at 150–250 μL per 20 mg tissue. Increase volume if lysis is incomplete; decrease volume to obtain higher protein concentration.
- c. Homogenize with a glass homogenizer until fully lysed. Alternatively, freeze tissue, grind in liquid nitrogen, then add lysis buffer.
- d. After complete lysis, centrifuge at 10,000–14,000 $\times g$ for 3–5 min. Collect the supernatant for downstream assays.

1. Lysis buffer volume guide: 150 μ L is usually sufficient per well of a 6-well plate or per 1 mL bacterial/yeast culture. For very high cell density, increase to 200 or 250 μ L. Lysate from 1×10^6 animal cells in 100 μ L buffer typically yields 2–4 mg/mL protein (varies by cell type).
2. Lysate from 20 mg frozen mouse liver in 200 μ L buffer typically yields 15–25 mg/mL protein (varies by tissue and condition).
3. Very fine tissue samples may be directly lysed after minimal chopping, using vigorous vortexing. This method is convenient but less thorough than homogenization or grinding.
4. For optimal performance, aliquot before use to avoid repeated freeze–thaw cycles.
5. All lysis steps must be performed on ice or at 4°C.
6. This product is for scientific research use only.



APEX^{BIO} Technology
www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com

