

2X Taq PCR Master Mix Plus (with dye)

Product description

Taq DNA polymerase can synthesize DNA according to the template under appropriate conditions, in the presence of gene-specific primers and dNTPs. Taq DNA polymerase has the ability of 5'→3' DNA polymerase and a weak 5'→3' exonuclease, but does not have the ability of 3'→5' exonuclease, which means that a dA protrusion will appear at the 3' end. Using this feature, it can be used for TA cloning.

This product contains the modified Taq DNA Polymerase, dNTP and an optimized buffer system. The Taq DNA Polymerase used has strong amplification performance and can efficiently amplify DNA fragments of ≤7 kb. It can achieve efficient amplification of different GC contents (30% - 70%). It is compatible with a variety of complex templates and is suitable for a wide range of applications, including plant and colony PCR as well as mouse gene identification. This Master Mix operation only requires the addition of primers, templates and ddH₂O, which greatly simplifies the experimental steps, reduces personal errors and improves the repeatability of the results. This product contains blue dye, which means it can be directly subjected to electrophoresis after amplification without the need to add sample loading buffer. If agarose Gel electrophoresis of PCR products is required, our product SYBR Safe DNA Gel Stain (No. A8743) can be added to the gel.

Composition and storage conditions

Components	Size						Storage
	1 mL	5×1 mL	20×1 mL	50×1 mL	100×1 mL		
2X Taq PCR Master Mix Plus (with dye)	1 mL	5×1 mL	20×1 mL	50×1 mL	100×1 mL		-20°C
Shipping: Dry Ice		Shelf life: 2 years					

Experimental operation

1. Reaction Setup

- Thaw all the components, mix them thoroughly and then centrifuge. Place on ice for later use.
- Configure the reaction system according to the following table:

Reagent	50 μ L Reaction	Final Concentration
2X Taq PCR Master Mix Plus (with dye)	25 μ L	1 \times
Forward primer (10 μ M)	1-2 μ L	0.2-0.4 μ M
Reverse primer (10 μ M)	1-2 μ L	0.2-0.4 μ M
Template DNA	X μ L	
ddH ₂ O	to 50 μ L	

***Note:**

- Generally, the final primer concentration is recommended to be 0.2-0.4 μ mol/L, but it can also be adjusted within the range of 0.1-1.0 μ mol/L as needed.
- The optimal reaction concentrations for different templates vary. Taking the 50 μ L system as an example, when the template is genomic DNA, the generally recommended template usage is 10-200 ng. When the template is plasmid or viral DNA, the generally recommended usage amount is 10 pg-5 ng. Excessive template quantity is prone to cause non-specific amplification.

- Invert and mix well, or gently pipette or gently vortex to thoroughly mix the reaction solution, and briefly centrifuge to collect the liquid to the bottom of the tube.

2. PCR Cycling Conditions

Place the reaction tube on the PCR instrument and start the thermal cycle. Set the reaction program according to the following table:

Procedure	Temperature	Time	Cycles
Initial Denaturation	95°C	3-5 min	1
Denaturation	95°C	30 s	30-35
Annealing	55-65°C	30 s	
Extension	72°C	15-30 s/kb	
Final extension	72°C	5 min	1

***Note:**

- When performing colony PCR, pre-denaturation for 10 min can fully break the cell walls, and both *Escherichia coli* and yeast can be efficiently amplified.
- The annealing temperature should be set according to the T_m value of the primer. If necessary, it is recommended to establish a temperature gradient to find the optimal temperature for the combination of primers and templates. In addition, the annealing temperature directly determines the amplification specificity. If poor amplification specificity is found, the annealing temperature can be appropriately increased.
- If the length of the target fragment is <3 kb, the extension time can be shortened to 15 s/kb. The length of the target fragment is > 3 kb,

and the recommended extension time is 30 s/kb. To achieve the best amplification effect or higher yield, it is recommended to uniformly extend at a speed of 30 s/kb.

Notes

Primer Design Guidelines: The primers for qPCR experimental design need to have good amplification effects and few non-specific products. The following design guidelines can be followed:

- 1) The last base at the 3' end of the primer is preferably G or C. The overall distribution of primers A, G, C and T should be as uniform as possible, and areas with high GC or AT content should be avoided.
- 2) The last 8 bases at the 3' end of the primer should avoid consecutive mismatches and the formation of hairpin structures.
- 3) It is preferable that the difference between the T_m values of the forward Primer and the reverse primer does not exceed 1°C, and the T_m value should be adjusted to 55-65° C (It is recommended to calculate the T_m value of the primer using Primer Premier 5).
- 4) Additional sequences of primers, that is, sequences that are not paired with the template, should not be involved in the calculation of primer T_m values. The GC content of the primers is controlled between 40% and 60%.
- 5) Avoid having complementary sequences of more than five bases within the primers or between the two primers, and avoid having complementary sequences of more than three bases at the 3' ends of the two primers.
- 6) After the primer design is completed, please use the NCBI BLAST function to search for primer specificity to avoid non-specific amplification.

APExBIO Technology

www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com