



5' DNA Adenylation Kit



Product description

5' DNA Adenylation Kit is specifically designed for the adenylation modification of 5'-phosphorylated single-stranded DNA (ssDNA) at the 5' end. This kit contains Mth RNA Ligase (an RNA ligase derived from thermophilic archaea) and ATP. During the adenylation reaction, 5'-phosphorylated single-stranded DNA (5'-pDNA) serves as the substrate. The Mth RNA Ligase, which possesses adenylation activity, hydrolyzes ATP into AMP and PPi, and subsequently transfers the AMP to the 5'-phosphate group of the ssDNA, ultimately generating 5'-adenylated ssDNA (5'-AppDNA).

This kit offers a simple and efficient one-step reaction for quantitative adenylation, typically converting 95% of pDNA to AppDNA, eliminating the need for product purification. It is more convenient than current chemical or enzymatic methods. The reaction temperature of 65°C effectively reduces secondary structure formation, and is suitable for reaction systems ranging from pmol to µmol substrate amounts. This product is commonly used to generate adenylated linkers for 3'-end ligation to ssRNA/ssDNA with a 3'-hydroxyl group, facilitating applications such as cloning, high-throughput sequencing library preparation, or PCR detection.

Composition and storage conditions

Size Size		Jie e Perectuii	
Components	10 rxns	50 rxns	Storage
Mth RNA Ligase	20 μL	100 μL	-20°C
10× Adenylation Reaction Buffer	100 μL	500 μL	-20°C
ATP (1 mM)	25 μL	125 μL	-20°C
Shinning: Dry Ice She	If life: 12 months		E Unitrovin

Experimental operation

Adenylation of the 5' end of phosphorylated ssDNA/ssRNA:

1. Refer to the table below to prepare the reaction system on ice:

Reagent	20 μL Reaction	Final Concentration
10× Adenylation Reaction Buffer	2 μL	1×

Phosphorylated ssDNA or ssRNA	ΧμΣ	5 μΜ
ATP (1 mM)	2 μL	0.1 mM
Mth RNA Ligase	2 μL	
Nuclease-free Water	to 20 μL	40

*Note:

- a. The adenylylation reaction system can be scaled up according to the proportions in the table above, with no significant effect on the yield of the adenylylated product.
- b. If RNA manipulation is involved, strict adherence to RNA handling protocols is required to avoid RNase contamination. Relevant reagents and consumables should be DEPC-treated to remove RNase or be confirmed as RNase-free. If ssRNA is involved, it is recommended to add an appropriate amount of RNase Inhibitor (Cat. No. K1046).
- 2. Incubate at 65°C for 60 min.

*Note: In some cases, to make the connection reaction more thoroughly, the reaction time can be appropriately extended.

- 3. Incubate at 85°C for 5 min to terminate the reaction.
- 4. Perform electrophoretic identification of adenylylated products. Before loading the samples for electrophoresis, incubate them at 95°C for 5 min and then place on ice to ensure complete denaturation. It is recommended to use a 15% denatured polyacrylamide gel containing 7 M urea for gel electrophoresis analysis.

*Note: The molecular weight of the sample will increase after adenosylation modification.

5. The sample can subsequently be concentrated using the conventional ethanol precipitation method, and then subjected to 3'-hydroxyl ligation with small RNAs using enzymes such as T4 RNA Ligase 2.

Notes

- 1. The substrate required for adenylylation must be ssDNA/ssRNA with a 5' phosphate group. The 3' end can be aminated or otherwise blocked, or it can remain unblocked.
- 2. The optimal reaction temperature for the enzyme included in the kit is 65°C. De-adenylylation may occur at 25°C. Therefore, it is recommended to heat the reaction at 85°C for 5 min after completion to inactivate the enzyme. If the enzyme is not inactivated, exposure to room temperature conditions in subsequent steps may lead to a decrease in the adenylylation ratio.
- 3. This product is for scientific use only!

Educe of Servery Light of the University

ARE Indicate the land of the l

ARE ENGLANDED

A PER Laboration Laboration International

ARE Language and Area of the A

ARE LEVE TO THE LE



APExBIO Technology

www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com