

Product Information

HAT Activity Fluorometric Assay Kit

I. Kit Contents:

Components	K2034-100 100 assays	Cap Color	Part Number
HAT Assay Buffer	25 ml	WM	K2034-C-1
Acetyl CoA (Lyophilized)	1 vial	Red	K2034-C-2
H3 Peptide (Lyophilized)	1 vial	Brown	K2034-C-3
Substrate Mix (Lyophilized)	1 vial	Green	K2034-C-4
Developer	100 µl	Orange	K2034-C-5
PicoProbe™	200 µl	Blue	K2034-C-6
CoA Standard (Lyophilized)	1 vial	Yellow	K2034-C-7
Positive Control (HeLa Nuclear Extract)	40 µl	Violet	K2034-C-8

II. Introduction:

Histone acetyltransferase (HAT) is an enzyme that acetylates conserved lysine on a histone protein by transferring an acetyl group from acetyl CoA. Histone acetylation can increase gene expression. HAT plays important roles in gene transcription, cell proliferation and differentiation.

The HAT Activity Fluorometric Assay Kit provides a sensitive and sensitive way for detection of HAT activity in a variety of samples based on fluorometric method that eliminates radioactivity in traditional assays. The assay utilizes H3 histone peptide and Acetyl CoA as substrates. The active HAT transfers acetyl groups from Acetyl-CoA to the H3 histone peptide, thereby producing two products - acetylated H3 histone peptide and CoA-SH. The CoA-SH reacts with the developer to yield a fluorophore product that can be easily detected using a fluorometer or a fluorescence plate reader at Ex/Em = 535/587 nm. The kit can determine HAT activity as low as 6 mU in a variety of samples.

III. Applications:

Measurement of HAT activity in Nuclear Extracts.

Measurement of HAT activity of purified enzyme preparations.

IV. Sample Type:

Nuclear extracts from cells and tissue.

Recombinant enzyme.

V. User Supplied Reagents and Equipment:

96-well plate with flat bottom. White plates are preferred for this assay.

Multi-well spectrophotometer capable of fluorescence detection.

VI. Storage and Handling:

Store kit at -80°C, protected from light. Warm Assay Buffer to room temperature before use. Briefly centrifuge small vials prior to opening. Read the entire protocol before performing the experiment.

VII. Reagent Preparation and Storage Conditions:

Acetyl CoA: Reconstitute with 220 μ l deionized water. Make 20 μ l aliquots and store at -80°C . Stable at -80°C for two months. Avoid repeated freeze/thaw. Keep on ice while in use.

H3 Peptide: Reconstitute with 420 μ l HAT Assay Buffer. Pipette up and down to dissolve completely. Aliquot and store at -80°C . Avoid repeated freeze/thaw. Use within two months. Keep on ice while in use.

Substrate Mix: Reconstitute with 1.1 ml HAT Assay Buffer. Pipette up and down to dissolve completely. Store at -80°C . Use within two months.

Developer: Store at -20°C . The solution is very viscous and difficult to pipette accurately. Immediately prior to use, take the required volume of developer and dilute 1:1 with an equal volume of HAT Assay Buffer.

PicoProbe™: Warm to room temperature and mix well before use. Store at -20°C .

CoA Standard: Reconstitute with 100 μ l HAT Assay Buffer to generate 100 mM solution & mix completely. Aliquot and store at -80°C . Avoid repeated freeze/thaw. Use within two months.

Positive Control: Aliquot & store at -80°C . Avoid repeated freeze/thaw. Use within two months.

VIII. HAT Activity Assay Protocol:

1. Sample Preparation: Prepare nuclear extract using BioVision's Nuclear/Cytosol Fractionation Kit. Add 2-10 μ l of sample and make up the volume to 50 μ l with HAT Assay Buffer. Add 50 μ l HAT Assay Buffer to one of the wells as Background Control. For Positive Control, add 2-4 μ l of HeLa Nuclear Extract into desired wells and make up the volume to 50 μ l with HAT Assay Buffer.

Note:

- For unknown samples, use varying sample amounts so as to obtain linear enzyme activity in the range of the Standard Curve.
- Dithiothreitol (DTT) and β -mercaptoethanol will interfere with the assay. Make sure samples are free of dithiothreitol (DTT) or β -mercaptoethanol.

2. Standard Curve Preparation: Dilute CoA Standard to 1 mM by adding 10 μ l of 100 mM CoA Standard to 990 μ l of HAT Assay Buffer. Dilute further to 0.1 mM by adding 10 μ l of 1 mM CoA Standard to 90 μ l of HAT Assay Buffer. Add 0, 2, 4, 6, 8 and 10 μ l of 0.1 mM CoA Standard into a series of wells in a 96-well plate to generate 0, 200, 400, 600, 800 and 1000 pmol/well of CoA Standard. Adjust the volume to 50 μ l/well with HAT Assay Buffer.

Note: Diluted CoA Standard is unstable. Discard the diluted Standard after use.

3. Reaction Mix: Mix enough reagents for the number of assays to be performed. Add reagents in the order shown. For each well, prepare 50 μ l Mix containing:

HAT Assay Buffer	30 μ l
H3 Peptide	4 μ l
Substrate Mix	10 μ l
Developer	2 μ l
PicoProbe™	2 μ l
Acetyl CoA	2 μ l

Add 50 μ l of the reaction mix to each well containing the Samples, Background Control, Standards and Positive Control. Mix well.

4. Measurement: Read fluorescence (Ex/Em = 535/587 nm) in kinetic mode at 25°C for 30-60 min. Choose two time points (T1 & T2) in the linear range of the plot and obtain the corresponding RFU for Sample (RS1 and RS2) and sample background (RB1 and RB2).

5. Calculation: Subtract 0 Standard reading from all Standard readings. Note: The CoA Standards will show some drift. Extrapolate the curve for each Standard to the Y-axis to obtain the Y-intercept. Plot the Standard Curve using the corrected intercept values. Calculate the HAT Activity of the test sample $\Delta\text{RFU} = (\text{RS2} - \text{RS1}) - (\text{RB2} - \text{RB1})$. Apply the ΔRFU to the Standard Curve to get B pmol of CoA formed during the reaction time ($\Delta\text{T} = \text{T2} - \text{T1}$).

Sample HAT Activity = $B/(\Delta\text{T} \times V) \times D = \text{pmol}/\text{min}/\text{ml} = \mu\text{U}/\text{ml}$

Where: B = CoA amount from Standard Curve (pmol)

ΔT = Reaction time (min.)

V = Sample volume added into the reaction well (ml)

D = Dilution Factor

Sample HAT Activity can also be expressed in $\mu\text{U}/\mu\text{g}$ of protein.

Unit Definition: One unit of HAT activity is the amount of enzyme that will generate 1.0 μmol of CoA per min. at 25°C using kit assay conditions.




Figure: (a) Co A Standard Curve. (b) HAT Activity in HeLa Nuclear Extract. (c) Specific Activity of HeLa Nuclear Extract and purified recombinant pCAF. Assays were performed following the kit protocol.

For research use only! Not to be used in humans.

Our promise

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For more details, please visit <http://www.apexbt.com/> or contact our technical team.

Tel: +1-(832)696-8203

Fax: +1-832-641-3177

Email: sales@apexbt.com