

Product Information

Red Active Caspase-9 Staining Kit

I. Kit Contents:

Component	K2055-25	K2055-100	Part Number
	25 assays	100 assays	
Red-LEHD-FMK	25 μ l	100 μ l	K2055-C-1
Wash Buffer	50 ml	2 x 100 ml	K2055-C-2
Z-VAD-FMK	10 μ l	10 μ l	A1902

II. Introduction:

Caspases activation is crucial in apoptosis. Red Caspase-9 Staining Kit is an easy and sensitive way of detecting activated caspases in living cells. This assay uses the Caspase-9 inhibitor, LEHD-FMK, coupled to sulfo-rhodamine (Red-LEHD-FMK) as a marker. In apoptotic cells, the cell permeable and nontoxic Red-LEHD-FMK irreversibly binds to activated Caspase-9.

III. Caspase-9 Assay Protocol:

A. Staining Procedure:

1. Induce apoptosis in cells (1×10^6 /ml) by desired method. Concurrently incubate a control culture without induction. An additional negative control can be prepared by adding the caspase inhibitor Z-VAD-FMK at 1 μ l/ml to an induced culture to inhibit caspase activation.
2. Aliquot 300 μ l each of the induced and control cultures into eppendorf tubes.
3. Add 1 μ l of Red-LEHD-FMK into each tube and incubate for 0.5 - 1 hour at 37°C incubator with 5% CO₂.
4. Centrifuge cells at 3000 rpm for 5 minutes and remove supernatant.
5. Resuspend cells in 0.5 ml of Wash Buffer, and centrifuge again.
6. Repeat Step 5.

Proceed to B, C, or D depending on methods of analysis.

B. Quantification by Flow Cytometry:

For flow cytometric analysis, resuspend cells in 300 μ l of Wash buffer. Put samples on ice. Analyzing samples by flow cytometry using the FL-2 channel.

C. Detection by Fluorescence Microscopy:

For fluorescence microscopic analysis, resuspend cells in 100 μ l Wash buffer. Put one drop of the cell suspension onto a microslide and cover with a coverslip. Observe cells under a fluorescence microscope using rhodamine filter. Caspase-9 positive cells appear to have brighter red signals, whereas Caspase-9 negative control cells show much weaker signal.

D. Analysis by Fluorescence Plate Reader:

For analysis with fluorescence plate reader, resuspend cells in 100 μ l Wash Buffer and then transfer the cell suspension to each well of the black microtiter plate. Measure the fluorescence intensity at Ex/Em = 540/570 nm (Note: Ex/Em=488/570 nm will also work, although it's not an optimal wavelength). For control, use wells containing unlabeled cells.

General Troubleshooting Guide:

Problems	Cause	Solution
High Background	<ul style="list-style-type: none"> • Cell density is higher than recommended • Cells were not washed well with wash buffer after staining • Cells were Incubated for extended period of time • Use of extremely confluent cells • Contaminated cells 	<ul style="list-style-type: none"> • Refer to data sheet and use the suggested cell number • Use the wash buffer provided, and as instructed in the datasheet • Refer to data sheets and incubate for exact times • Perform assay when cells are at 70-95% confluency • Check for bacteria/ yeast/ mycoplasma contamination
Lower signal levels	<ul style="list-style-type: none"> • Cells did not initiate apoptosis • Very few cells used for analysis • Incorrect setting of the equipment used to read samples • Use of expired kit or improperly stored reagents 	<ul style="list-style-type: none"> • Determine the time-point for initiation of apoptosis after induction (time-course experiment) • Refer to data sheet for appropriate cell number • Refer to data sheet and use the recommended filter setting • Always check the expiry date and store the components appropriately
Erratic results	<ul style="list-style-type: none"> • Old (unhealthy) cells used • Adherent cells were dislodged and washed away prior to assaying • Incorrect incubation times or temperatures • Incorrect volumes used 	<ul style="list-style-type: none"> • Seed healthy cells and make sure cells are healthy prior to induction of apoptosis • Collect all cells (both attached and dislodged) after induction for accurate results • Refer to datasheet & verify correct incubation times and temperatures • Use calibrated pipettes and aliquot correctly
<p>Note: The most probable cause is listed under each section. Causes may overlap with other sections</p>		

For research use only! Not to be used in humans.

Our promise

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For more details, please visit <http://www.apexbt.com/> or contact our technical team.

Tel: +1-(832)696-8203

Fax: +1-832-641-3177

Email: sales@apexbt.com