

## Product Information

### ATP Cell Viability Bioluminescence Assay Kit

#### I. Kit Contents:

Components	K2102-200 200 assays	Cap Color	Part Number
Nucleotide Releasing Buffer	50 ml	NM	K2102-C-1
ATP Monitoring Enzyme	Lyophilized	Green	K2102-C-2
ADP Converting Enzyme	Lyophilized	Blue	K2102-C-3
Enzyme Reconstitution Buffer	2 ml	Red	K2102-C-4

#### II. Introduction:

The changes in ADP/ATP ratio is used to recognized the different states of the cell death and viability. Levels of ATP and ADP levels are increased in proliferating cells but decreased in apoptotic cells. The ADP/ATP Ratio Assay Kit detects ATP and ADP levels by bioluminescent for a fast screening of apoptosis, necrosis, growth arrest, and cell proliferation simultaneously in mammalian cells. This kit utilizes luciferase to catalyze the formation of light from ATP and luciferin, and the light can measured using a luminometer or Beta Counter. ADP level is measured by its conversion to ATP that is subsequently detected using the same reaction. The assay can be fully automatic for high throughput and is highly sensitive (detects 100 mammalian cells/well).

#### III. Reagent Reconstitution and General Consideration:

1. Reconstitute ATP Monitoring Enzyme with 2.1 ml of the Enzyme Reconstitution Buffer & mix gently by inversion.
2. Reconstitute ADP Converting Enzyme with 220  $\mu$ l of the Nucleotide Releasing Buffer & mix gently by inversion.
3. The reconstituted enzymes are stable for up to 2 months at 4 $^{\circ}$ C after reconstitution. For more accurate handling, the reconstituted ADP-Converting Enzyme can be diluted 10-fold with Nucleotide Releasing Buffer just before use (Section IV Step 7), then use 10  $\mu$ l of the enzyme for each assay.
4. The kit is significantly more sensitive than other methods used for cell viability assays. As a rule, we recommend using  $1 \times 10^4$  -  $1 \times 10^6$  cells per assay. Avoid contamination with ATP from exogeneous biological sources, such as bacteria or fingerprints.
5. Ensure that the Nucleotide Releasing Buffer is at room temperature before use. The optimal temperature is 22 $^{\circ}$ C. Keep enzymes on ice during the assay and protect from light as much as possible.
6. The assay can be performed using either a single tube or a white walled 96-well luminometer plate (100  $\mu$ l/well culture volume is recommended).

#### IV. ADP/ATP Ratio Assay Protocol:

1. Induce apoptosis in cells by desired method. Concurrently incubate a control culture without induction.
2. For each sample well to be measured, mix 100  $\mu$ l of reaction mix consisting of: ATP monitoring enzyme: 10  $\mu$ l Nucleotide Releasing Buffer: 90  $\mu$ l
3. Add 100  $\mu$ l of the reaction mix to the appropriate wells of a 96-well plate and read the background luminescence (Data A). For higher accuracy let the reaction mix sit at room temperature to burn off low level ATP contamination for a few hours.
4. For suspension cells, transfer 10  $\mu$ l of the cultured cells ( $10^3$  –  $10^4$ ) into luminometer plate.
5. For adherent cells, remove culture medium and treat cells ( $10^3$  –  $10^4$ ) with 50  $\mu$ l of Nucleotide Releasing Buffer for 5 minutes at room temperature with gentle shaking. Transfer into luminometer plate.
6. After ~ 2 minutes, read the sample in a luminometer or luminescence-capable plate reader (Data B).

7. To measure ADP levels in the cells, read the samples (from step 6) again (Data C), then add 1 µl of ADP Converting Enzyme. Read the samples after ~ 2 minutes (Data D).

Note: The results can be analyzed using cuvette-based luminometers or Beta Counters. When Beta Counter is used it should be programmed in the “out of coincidence” (or Luminescence mode) for measurement. The entire assay can directly be done in a 96-well plate. It can also be programmed automatically using instrumentation with injectors (When using injector the ATP Monitoring Enzyme and the ADP Converting Enzyme can be diluted with the Nuclear Releasing Buffer at 1:50 for injector).

The assay utilizes a “glow-type” luciferase which has replaced the original “flash-type” luciferase. While still sensitive to sub-picomole amounts of ATP, the glow-type reactions can still be read an hour later. This means that ATP & ADP levels are now determined by quasi-steady-state light output levels. This makes the reading of an entire 96-well (384-well) plate much more feasible.

8. Calculation:

ADP/ATP Ratio is calculated as:  $(\text{Data D} - \text{Data C}) / (\text{Data B} - \text{Data A})$

Interpretation of Results:

Cell Fate	ADP Level	ATP Level	ADP/ATP
Proliferation	Very Low	High	Very Low
Growth Arrest	Low	Slightly Increased	Low
Apoptosis	High	Low	High
Necrosis	Much Higher	Very Low	Much Higher

The interpretation of different ratios obtained may vary significantly according to the cell types and conditions used. However, the following criteria may be used as guidelines:

- Test gives markedly elevated ATP values with no significant increase in ADP levels in comparison to control cells = proliferation.
- Test gives similar or slightly higher levels of ATP and with little or no change in ADP compared to control = growth arrest.
- Test gives lower levels of ATP to control but shows an increase in ADP = apoptosis.
- Test gives considerable lower ATP levels than control but greatly increased ADP = necrosis.

**For research use only! Not to be used in humans.**

## Our promise

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For more details, please visit <http://www.apexbt.com/> or contact our technical team.

Tel: +1-(832)696-8203

Fax: +1-832-641-3177

Email: [sales@apexbt.com](mailto:sales@apexbt.com)