

Product Information

Glutathione Colorimetric Assay Kit

I. Kit Contents:

Components	K2106-100 100 assays	Cap Color	Part Number
Glutathione Reaction Buffer	100 ml	NM	K2106-C-1
Glutathione Substrate (DTNB)	2 vials	Red	K2106-C-2
NADPH Generating Mix (lyophilized)	2 vials	Blue	K2106-C-3
Glutathione Reductase (lyophilized)	2 vials	Green	K2106-C-4
Sulfosalicylic Acid (SSA, 1 gram)	1 bottle	WM	K2106-C-5
GSH Standard (lyophilized, MW 307)	2 x 1 mg	Yellow	K2106-C-6

II. Introduction:

Glutathione (GSH) is a low-molecular-weight thiol that play a vital role in protecting cell from oxidative stress in mammalian cells. The Glutathione Colorimetric Assay Kit provides an easy and colorimetric method for analyzing total glutathione or the reduced form of glutathione by a microtiter plate reader. The assay utilizes glutathione recycling system by DTNB and glutathione reductase. The yellow color 2-nitro-5-thiobenzoic acid that generate by the reaction of DTNB and GSH. Thus the GSH concentration can be measured at 412 nm absorbance. The generated GSSG can be reduced back to GSH by glutathione reductase, and GSH reacts with DTNB again to produce more 2-nitro-5-thiobenzoic acid. As a result, the recycling system significantly improves the sensitivity of total GSH detection. The 5-Sulfosalicylic acid (SSA) in the kit removes proteins from samples to protect the GSH oxidation and γ -glutamyl transpeptidase reaction. The kit can quantify glutathione from 1-100 ng/well in a 200 μ l reaction. Increase the reaction time to generate stronger signal for detecting lower glutathione concentrations, such as in blood samples. The kit can also specifically detect the reduced form of glutathione (GSH) by omitting the glutathione reductase from the reaction mixture. The sensitivity for detecting the reduced form of glutathione (without recycling system) is 100 times lower than detecting the total glutathione.

III. Sample Preparation:

Note: Peptide thiol may interfere with the assay of reduced form glutathione. SSA treatment may not able to complete remove all small peptides from samples. Further purification may be required to accurately measure reduced form glutathione. Peptide thiols don't significantly interfere with total glutathione assay.

A. Cell Sample Preparation (0.5 - 1 x 10⁶ cells/assay)

1. Treat cells by desired method. Concurrently incubate a control culture without treatment.
2. Collect cells by centrifugation at 700 x g for 5 minutes at 4°C Remove supernatant.
3. Resuspend cell pellet in 0.5 ml ice-cold PBS. Transfer into a 1.5 ml microcentrifuge tube, and centrifuge at 700 x g for 5 minutes at 4°C Remove supernatant.
4. Lyse cells in 80 μ l ice-cold Glutathione Buffer. Incubate on ice for 10 minutes.
5. Add 20 μ l of 5% SSA (see below for SSA preparation), mix well and centrifuge at 8000 x g for 10 min. Transfer supernatant to a fresh tube and use it for glutathione assay.

B. Tissue Sample Preparation (100 mg)

1. Homogenize the tissue in 0.4 ml of Glutathione Buffer.
2. Add 100 μ l of 5% SSA (see below for SSA preparation), mix well, and centrifuge at 8000 x g for 10 minutes.

3. Transfer supernatant to a fresh tube and use it for glutathione assay.

C. Plasma Sample Preparation

1. Centrifuge an anticoagulant treated blood at 1000 x g for 10 min at 4°C.
2. Transfer the top plasma layer to a new tube and add 1/4 vol of 5 % SSA. Mix well.
3. Centrifuge at 8000 x g for 10 min at 4°C.
4. Transfer supernatant to a new tube, and use it for the glutathione assay.

D. Erythrocyte Sample Preparation

1. Centrifuge an anticoagulant treated blood at 1000 x g for 10 min at 4°C.
2. Discard the supernatant and the white buffy layer.
3. Lyse the erythrocytes with 4 vol of Glutathione Buffer. Keep on ice for 10 min.
4. Add 1 vol 5% SSA, mix well, and centrifuge at 8000 x g for 10 minutes. Transfer supernatant to a fresh tube and use it for glutathione assay.

Note: Erythrocytes can be isolated from the remaining sample solution after the plasma sample isolation.

IV. Preparation of Solutions & Storage Conditions:

Substrate: Add 1 ml of Glutathione Buffer to 1 vial of substrate and dissolve it completely. Store the remaining solution at -20°C stable for 2 months.

NADPH Generating Mix: Add 1 ml of Glutathione Buffer to 1 vial of the NADPH mix. Store the solution at -20°C stable for 2 months.

Glutathione Reductase: Add 1 ml of Glutathione Buffer to 1 vial of the enzyme and dissolve. Use up the solution within 1 day. SSA: Add 19 ml of dH₂O to make 5% solution and then dilute 5 ml of the solution with Glutathione Buffer to make 1% SSA solution. Store at 4°C stable for 6 months.

GSH Standard: Add 1 ml of 1% SSA to the GSH standard vial to generate 1µg/µl GSH standard solution. Store at -20°C stable for 2 months.

V. Preparation of Solutions for Standard Curve:

To generate standard curve for detecting the reduced form of glutathione only, add 50, 40, 30, 20, 10, and 0 µl of the 1 µg/µl GSH standard into each labeled microcentrifuge tubes, add 1% SSA to make up for a total volume of 100 µl/tube.

To generate standard curve for detecting the total glutathione, dilute the 1 µg/µl glutathione solution into 10 ng/µl with 1% SSA. Add 50, 40, 30, 20, 10, and 0 µl of the 10 ng/µl GSH standard into each labeled microcentrifuge tubes, add 1% SSA to make up for a total volume of 100 µl/tube.

VI. Glutathione Assay Protocol:

1. Prepare enough Reaction Mix for the standard and samples to be assayed in 96-well plate (not provided). Each well should contain:

NADPH Generating Mix	20 µl
Glutathione Reductase	20 µl
Glutathione Reaction Buffer	120 µl

For detecting the reduced form of glutathione only, omit Glutathione Reductase. Use 20 µl of the Glutathione Reaction Buffer replace the 20 µl of glutathione Reductase.

2. Mix well. Add 160 µl of the Reaction Mix to each well and incubate at room temperature for 10 minutes to generate NADPH.

3. Add 20 µl of either the GSH standard solutions or the sample solution. Incubate the plate at room temperature for 5 - 10 min.

Note: We recommend to make several dilutions of your sample using the 1% SSA to make sure the readings are within the range of the standard calibration curve.

4. Add 20 µl of Substrate solution, and incubate at room temperature for 5-10 min, or longer if the samples contain low levels of glutathione.

Notes: a) Since the reaction starts immediately after the addition of substrate, use a multichannel pipette or repeating pipette is recommended to avoid the reaction time lag among wells.

b) You can read samples immediately and at various times following addition of the substrate solution for kinetic studies.

5. Read the absorbance at 405 nm or 415 nm using a microplate reader.

6. Determine concentrations of GSH in the sample solutions using the standard glutathione calibration curve.

Note:

A. Using reduced form glutathione Standard Curve for detecting reduced form of glutathione. Using total Glutathione Standard Curve for detecting total glutathione. There are about 10 to 100 fold difference in detection sensitivity between detecting reduced form glutathione and total glutathione (see procedure step IV for preparation of standard curve).

B. The colorimetric reaction is stable and the O.D. increases linearly over 30 min for total glutathione detection.

VII. Calculation of Total Glutathione:

Pseudo-end point method: Total Glutathione = $(O.D._{sample} - O.D._{blank})$

Kinetic method: Total Glutathione = $(Slope_{sample} - Slope_{blank}) / Slope$.

VIII. Reagent Interference:

Reducing agents such as ascorbic acid, β -mercaptoethanol, dithiothreitol (DTT) and cysteine, or thiol reactive compounds such as maleimide compounds, interfere with the glutathione assay and therefore should be avoided during the sample preparation. When detecting the reduced form of glutathione, protein thiols can generate significant background signal. In such cases, it is necessary to completely remove proteins from samples. We suggest using Centrifugal Spin column with 10 kDa molecular weight cut off filter to remove proteins. Then the reduced glutathione can be easily detected from spin through samples.




Fig .2. Glutathione Standard Curve. Various amounts of standard glutathione was added to the glutathione reaction and incubated for 10 min according to the kit instructions. Absorbance was measured at O.D_{405 nm}.

General Troubleshooting Guide:

Problems	Cause	Solution
Assay not working	<ul style="list-style-type: none"> • Use of ice-cold reaction buffer • Omission of a step in the protocol • Plate read at incorrect wavelength • Use of a different 96-well plate 	<ul style="list-style-type: none"> • Reaction buffer must be at room temperature • Refer and follow the data sheet precisely • Check the wavelength in the data sheet and the filter settings of the instrument • Fluorescence: Black plates ; Luminescence: White plates; Colorimeters: Clear plates

Samples with erratic readings	<ul style="list-style-type: none"> • Use of an incompatible sample type • Samples prepared in a different buffer • Samples were not deproteinized (if indicated in datasheet) • Cell/ tissue samples were not completely homogenized • Samples used after multiple free-thaw cycles • Presence of interfering substance in the sample • Use of old or inappropriately stored samples 	<ul style="list-style-type: none"> • Refer data sheet for details about incompatible samples • Use the reaction buffer provided in the kit or refer data sheet for instructions • Use the 10 kDa spin cut-off filter or PCA precipitation as indicated • Use Dounce homogenizer (increase the number of strokes); observe for lysis under microscope • Aliquot and freeze samples if needed to use multiple times • Troubleshoot if needed, deproteinize samples • Use fresh samples or store at correct temperatures till use
Lower/ Higher readings in Samples and Standards	<ul style="list-style-type: none"> • Improperly thawed components • Use of expired kit or improperly stored reagents • Allowing the reagents to sit for extended times on ice • Incorrect incubation times or temperatures • Incorrect volumes used 	<ul style="list-style-type: none"> • Thaw all components completely and mix gently before use • Always check the expiry date and store the components appropriately • Always thaw and prepare fresh reaction mix before use • Refer data sheet & verify correct incubation times and temperatures • Use calibrated pipettes and aliquot correctly
Readings do not follow a linear pattern for Standard curve	<ul style="list-style-type: none"> • Use of partially thawed components • Pipetting errors in the standard • Pipetting errors in the reaction mix • Air bubbles formed in well • Standard stock is at an incorrect concentration • Calculation errors • Substituting reagents from older kits/ lots 	<ul style="list-style-type: none"> • Thaw and re-suspend all components before preparing the reaction mix • Avoid pipetting small volumes • Prepare a master reaction mix whenever possible • Pipette gently against the wall of the tubes • Always refer the dilutions in the data sheet • Recheck calculations after referring the data sheet • Use fresh components from the same kit
Unanticipated results	<ul style="list-style-type: none"> • Measured at incorrect wavelength • Samples contain interfering substances • Use of incompatible sample type • Sample readings above/below the linear range 	<ul style="list-style-type: none"> • Check the equipment and the filter setting • Troubleshoot if it interferes with the kit • Refer data sheet to check if sample is compatible with the kit or optimization is needed • Concentrate/ Dilute sample so as to be in the linear range

Note: The most probable list of causes is under each problem section. Causes/ Solutions may overlap with other problems.

For research use only! Not to be used in humans.

Our promise

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For more details, please visit <http://www.apexbt.com/> or contact our technical team.



Tel: +1-(832)696-8203

Fax: +1-832-641-3177

Email: sales@apexbt.com