

Product Information

Caspase Fluorometric Substrate Set Plus

I. Kit Contents:

Description	Concentration	7 x 25 assays	Part Number
Caspase-1 Substrate, Ac-YVAD-AFC	1 mM	125 µl	A8956
Caspase-2 Substrate, Ac-VDVAD-AFC	1 mM	125 µl	K2146-C-1
Caspase-3 Substrate, Ac-DEVD-AFC	1 mM	125 µl	K2146-C-2
Caspase-5 Substrate, Ac-WEHD-AFC	1 mM	125 µl	K2146-C-3
Caspase-6 Substrate, Ac-VEID-AFC	1 mM	125 µl	K2146-C-4
Caspase-8 Substrate, Ac-IETD-AFC	1 mM	125 µl	K2146-C-5
Caspase-9 Substrate, Ac-LEHD-AFC	1 mM	125 µl	A4419
Cell Lysis Buffer	N/A	100 ml	K2146-C-6
2X Reaction Buffer	N/A	20 ml	K2146-C-7
DTT	1 M	0.4 ml	K2146-C-8

II. Introduction:

Caspases (Cysteine-dependent aspartate-directed proteases) are a family of cysteine proteases that play important roles in apoptosis, inflammation and necrosis. Sequential activation of caspases plays a critical role in cell apoptosis. Caspase-1, 2, 3, 4, 5, 6, 7, 8, 9 and 10 are members of caspases family. Caspase-1, 4 and 5 are involved in T-cell maturation. Caspase-2, 8, 9 and 10 are initiator caspases that cleave inactive pro-forms of effector caspases into active caspases. Caspase-3, 6 and 7 are effector caspases that cleave other protein substrates within the cell to trigger the apoptotic process.

Caspase Fluorometric Substrate Set Plus is composed of ready-to-use AFC-labeled substrates for caspase-1, 2, 3, 5, 6, 8 and 9 of caspase family proteases. The kit is used to detect activities of members of caspase family proteases. All substrates are provided in liquid ready-to-use form.

III. Assay Procedure:

1. Induce apoptosis in cells by desired method. Concurrently incubate a control culture without induction.
2. Count cells and pellet $1 - 5 \times 10^6$ cells or use 50 - 200 µg cell lysates if protein concentration has been measured.
3. Resuspend cells in 50 µl of chilled Cell Lysis Buffer.
4. Incubate cells on ice for 10 minutes.
5. Add 50 µl of 2X Reaction Buffer and 1 µl DTT to each sample.
6. Add 5 µl of the 1 mM AFC conjugated substrates (50 µM final conc.) into each tube individually and incubate at 37°C for 1 - 2 hour.
7. Read samples in a fluorometer equipped with a 400 nm excitation filter and 505 nm emission filter. For a plate-reading set-up, transfer the samples to a 96-well plate. You may perform the entire assay directly in a 96-well plate.

Fold-increase in caspase activity can be determined by comparing these results with the level of the uninduced control.

General Troubleshooting Guide for Caspase Colorimetric and Fluorometric Kits:

Problems	Cause	Solution
Assay not working	<ul style="list-style-type: none"> • Cells did not lyse completely • Experiment was not performed at optimal time after 	<ul style="list-style-type: none"> • Resuspend the cell pellet in the lysis buffer and incubate as described in the datasheet

	<p>apoptosis induction</p> <ul style="list-style-type: none"> • Plate read at incorrect wavelength • Old DTT used 	<ul style="list-style-type: none"> • Perform a time-course induction experiment for apoptosis • Check the wavelength listed in the datasheet and the filter settings of the instrument • Always use freshly thawed DTT in the cell lysis buffer
High Background	<ul style="list-style-type: none"> • Increased amount of cell lysate used • Increased amounts of components added due to incorrect pipetting • Incubation of cell samples for extended periods • Use of expired kit or improperly stored reagents • Contaminated cells 	<ul style="list-style-type: none"> • Refer to datasheet and use the suggested cell number to prepare lysates • Use calibrated pipettes • Refer to datasheet and incubate for exact times • Always check the expiry date and store the individual components appropriately • Check for bacteria/ yeast/ mycoplasma contamination
Lower signal levels	<ul style="list-style-type: none"> • Cells did not initiate apoptosis • Very few cells used for analysis • Use of samples stored for a long time • Incorrect setting of the equipment used to read samples • Allowing the reagents to sit for extended times on ice 	<ul style="list-style-type: none"> • Determine the time-point for initiation of apoptosis after induction (time-course experiment) • Refer to datasheet for appropriate cell number • Use fresh samples or aliquot and store and use within one month for the assay • Refer to datasheet and use the recommended filter setting • Always thaw and prepare fresh reaction mix before use
Samples with erratic readings	<ul style="list-style-type: none"> • Uneven number of cells seeded in the wells • Samples prepared in a different buffer • Adherent cells dislodged and lost at the time of experiment • Cell/ tissue samples were not completely homogenized • Samples used after multiple freeze-thaw cycles • Presence of interfering substance in the sample • Use of old or inappropriately stored samples 	<ul style="list-style-type: none"> • Seed only equal number of healthy cells (correct passage number) • Use the cell lysis buffer provided in the kit • Perform experiment gently and in duplicates/triplicates; apoptotic cells may become floaters • Use Dounce homogenizer (increase the number of strokes); observe efficiency of lysis under microscope • Aliquot and freeze samples, if needed to use multiple times • Troubleshoot as needed • Use fresh samples or store at correct temperatures until use
Unanticipated results	<ul style="list-style-type: none"> • Measured at incorrect wavelength • Cell samples contain interfering substances 	<ul style="list-style-type: none"> • Check the equipment and the filter setting • Troubleshoot if it interferes with the kit (run proper controls)
General issues	<ul style="list-style-type: none"> • Improperly thawed components • Incorrect incubation times or temperatures • Incorrect volumes used • Air bubbles formed in the well/tube • Substituting reagents from older kits/ lots • Use of a different 96-well plate 	<ul style="list-style-type: none"> • Thaw all components completely and mix gently before use • Refer to datasheet & verify the correct incubation times and temperatures • Use calibrated pipettes and aliquot correctly • Pipette gently against the wall of the well/tubes • Use fresh components from the same kit • Fluorescence: Black plates; Absorbance: Clear plates

Note# The most probable cause is listed under each section. Causes may overlap with other sections.



For research use only! Not to be used in humans.

Our promise

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For more details, please visit <http://www.apexbt.com/> or contact our technical team.

Tel: +1-(832)696-8203

Fax: +1-832-641-3177

Email: sales@apexbt.com