

Product Information

Caspase Colorimetric Substrate Set II Plus

I. Kit Contents:

Description	Concentration	25 assays	
Caspase-1 Substrate, Ac-YVAD-pNA	4 mM	125 µl	
Caspase-2 Substrate, Ac-VDVAD-pNA	4 mM	125 µl	
Caspase-3 Substrate, Ac-DEVD-pNA	4 mM	125 µl	
Caspase-4 Substrate, Ac-LEVD-pNA	4 mM	125 µl	
Caspase-5 Substrate, Ac-WEHD-pNA	4 mM	125 µl	
Caspase-6 Substrate, Ac-VEID-pNA	4 mM	125 µl	
Caspase-8 Substrate, Ac-IETD-pNA	4 mM	125 µl	
Caspase-9 Substrate, Ac-LEHD-pNA	4 mM	125 µl	
Caspase-10 Substrate, Ac-AEVD-pNA	4 mM	125 µl	
Cell Lysis Buffer	N/A	100 ml	
Dilution Buffer	N/A	200 ml	
2X Reaction Buffer	N/A	20 ml	
DTT	1 M	0.4 ml	

II. Introduction:

Caspases (Cysteine-dependent aspartate-directed proteases) are a family of cysteine proteases that play important roles in apoptosis, inflammation and necrosis. Sequential activation of caspases plays a critical role in cell apoptosis. Caspase-1, 2, 3, 4, 5, 6, 7, 8, 9 and 10 are members of caspases family. Caspase-1, 4 and 5 are involved in T-cell maturation. Caspase-2, 8, 9 and 10 are initiator caspases that cleave inactive pro-forms of effector caspases into active caspases. Caspase-3, 6 and 7 are effector caspases that cleave other protein substrates within the cell to trigger the apoptotic process.

Caspase Colorimetric Substrate Set II Plus is composed of ready-to-use pNA-labeled substrates for Caspase-1, 2, 3, 4, 5, 6, 8, 9 and 10 of caspase family proteases. The kit is used to detect activities of members of caspase family proteases. All substrates are provided in liquid form.

III. Assay Procedure:

1. Induce apoptosis in cells by desired method. Concurrently incubate a control culture without induction.

2. Count cells and pellet 1 - 5 x 106 cells.

3. Resuspend cells in 50 μ l of chilled Cell Lysis Buffer and incubate cells on ice for 10 minutes. Centrifuge for 1 min in a microcentrifuge (10,000 x g).

- 5. Transfer supernatant to a fresh tube and assay protein Concentration.
- 6. Dilute 100 300 μg protein to 50 μl Cell Lysis Buffer for each assay.
- 7. Add 50 µl of 2X Reaction Buffer containing 10 mM DTT to each sample.
- 8. Add 5 µl of the 4 mM pNA conjugated substrates (200 µM final conc.) into each tube individually and incubate at 37 °C for 1 2 hour.

9 Read samples at 400 or 405 nm in a microtiter plate reader, or spectrophotometer using a 100 μ l micro quartz cuvette, or dilute sample to 1 ml with Dilution Buffer and using regular cuvette (note: Dilution of the samples proportionally decreases the reading).

Fold-increase in caspase activity can be determined by comparing these results with the level of the uninduced control.



Note: Background reading from cell lysates and buffers must be subtracted from the readings of both induced and the uninduced samples before you calculate the fold increase in caspase activity.

Problems	Cause	Solution
Assay not working	• Cells did not lyse completely	• Resuspend the cell pellet in the lysis buffer and incubate as
	• Experiment was not performed at optimal time after	described in the datasheet
	apoptosis induction	• Perform a time-course induction experiment for apoptosis
	Plate read at incorrect wavelength	• Check the wavelength listed in the datasheet and the filter
	• Old DTT used	settings of the instrument
		• Always use freshly thawed DTT in the cell lysis buffer
High Background	Increased amount of cell lysate used	• Refer to datasheet and use the suggested cell number to
	• Increased amounts of components added due to incorrect	prepare lysates
	pipetting	• Use calibrated pipettes
	• Incubation of cell samples for extended periods	• Refer to datasheet and incubate for exact times
	• Use of expired kit or improperly stored reagents	· Always check the expiry date and store the individual
	Contaminated cells	components appropriately
		Check for bacteria/ yeast/ mycoplasma contamination
Lower signal	Cells did not initiate apoptosis	• Determine the time-point for initiation of apoptosis after
levels	• Very few cells used for analysis	induction (time-course experiment)
	• Use of samples stored for a long time	• Refer to datasheet for appropriate cell number
	• Incorrect setting of the equipment used to read samples	• Use fresh samples or aliquot and store and use within one
	• Allowing the reagents to sit for extended times on ice	month for the assay
		• Refer to datasheet and use the recommended filter setting
		• Always thaw and prepare fresh reaction mix before use
Samples with	• Uneven number of cells seeded in the wells	• Seed only equal number of healthy cells (correct passage
erratic readings	• Samples prepared in a different buffer	number)
	• Adherent cells dislodged and lost at the time of experiment	• Use the cell lysis buffer provided in the kit
	Cell/ tissue samples were not completely homogenized	• Perform experiment gently and in duplicates/triplicates;
	Samples used after multiple freeze-thaw cycles	apoptotic cells may become floaters
	• Presence of interfering substance in the sample	• Use Dounce homogenizer (increase the number of strokes);
	• Use of old or inappropriately stored samples	observe efficiency of lysis under microscope
		• Aliquot and freeze samples, if needed to use multiple times
		Troubleshoot as needed
		• Use fresh samples or store at correct temperatures until use
Unanticipated	Measured at incorrect wavelength	• Check the equipment and the filter setting
results	Cell samples contain interfering substances	• Troubleshoot if it interferes with the kit (run proper
		controls)
General issues	Improperly thawed components	• Thaw all components completely and mix gently before use
	Incorrect incubation times or temperatures	• Refer to datasheet & verify the correct incubation times and
	Incorrect volumes used	temperatures
	• Air bubbles formed in the well/tube	• Use calibrated pipettes and aliquot correctly

General	Troubles	shooting	Guide for	Caspase	Colorimetric	and Fluorom	etric Kits:



	• Substituting reagents from older kits/ lots	• Pipette gently against the wall of the well/tubes	
	• Use of a different 96-well plate	• Use fresh components from the same kit	
		• Fluorescence: Black plates; Absorbance: Clear plates	
Note# The most probable cause is listed under each section. Causes may overlap with other sections.			

For research use only! Not to be used in humans.

Our promise

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For more details, please visit <u>http://www.apexbt.com/</u> or contact our technical team.

Tel: +1-(832)696-8203 Fax: +1-832-641-3177 Email: sales@apexbt.com