

Caspase-5 Fluorometric Assay Kit

Introduction

Cysteine-containing aspartate proteolytic enzymes (Caspase) are a family of cysteine proteases that play important roles in apoptosis, necrosis, and inflammation. Caspase-1/4/5/11 are very important in inflammatory response and pyroptosis. Caspase-4/5 can recognize and bind intracellular lipopolysaccharide (LPS), which subsequently activate inflammasomes, which in turn promote the maturation and release of inflammatory mediator IL-1 β , leading to an inflammatory response. In addition, activated Caspase-4/5 can also cleave the substrate protein GSDMD, triggering the formation of cell membrane pores, resulting in the release of cell contents, and ultimately pyroptosis.

Caspase-5 Fluorometric Assay Kit provides a convenient and simple way for detecting the WEHD-dependent caspase activity. WEHD-AFC (AFC:7-amino-4-trifluoromethyl coumarin) emits blue light ($\lambda_{\text{max}} = 400 \text{ nm}$); while cleavage of WEHD-AFC by Caspase-5 or related caspases, free AFC emits a yellow-green fluorescence ($\lambda_{\text{max}} = 505 \text{ nm}$), which can be quantified by using a microplate reader or a fluorometer. Comparison of the fluorescence of AFC from an apoptotic sample with an uninduced control determines the fold increase in Caspase-5 activity.

Components and Storage

Components	K2195-25T	K2195-100T	K2195-200T	K2195-400T
Cell Lysis Buffer	25 mL	100 mL	100 mL	100 mL
2X Reaction Buffer	2 mL	4 X 2 mL	16 mL	32 mL
WEHD-AFC (1 mM)	125 μL	500 μL	2X 0.5 mL	2X 1 mL
DTT (1 M)	100 μL	400 μL	400 μL	400 μL

Store the kit at -20°C , stable for 6 months. WEHD-AFC (1 mM) should be stored away from light.

Protocol

1. Preparation before the experiment:

- 1) Equilibrate Cell Lysis Buffer and 2X Reaction Buffer to room temperature before use. Thaw WEHD-AFC (1 mM) in ice before use. Once thawed, aliquot WEHD-AFC (1 mM) to avoid repeated freeze/thaw cycles.
- 2) Add 10 μL DTT (1 M) per 1 mL of 2X Reaction Buffer to make the 2X Reaction Buffer (containing DTT).

2. Sample Preparation:

Try to use fresh samples. If this is not possible in time, it is recommended to complete the sample preparation step before storing at -80°C . And thaw and mix them on ice before use.

1) For cell samples:

- a) Prepare two groups of cells. One experimental group is treated according to the experimental design and the other is an untreated control group.
- b) Collect $1-5 \times 10^6$ cells after treatment. Wash the cells one time with pre-chilled PBS.

***Note:** When using adherent cells for testing, if there are suspended cells after treatment in the experimental group, the suspended cells should be collected and tested together.

- c) Resuspend cells in 50 μL of pre-chilled Cell Lysis Buffer. Incubate in ice for 10 min.
- d) After incubation, centrifuge samples at 12,000 rpm for 2-5 min, carefully transfer the supernatant to a new EP tube and place on ice for later use.
- e) Take a small amount of supernatant (1-2 μL) to measure the protein concentration Cpr using the Bradford method. To ensure the accuracy of the experiment, it is recommended to achieve a protein concentration of 1-4 mg/mL.

2) For tissue samples:

- a) Harvest 50-100 mg of tissue samples from the experimental group and the untreated control group. Wash cells with pre-chilled PBS once.
- b) Add 500 μL of pre-chilled Cell Lysis Buffer to homogenize tissue on ice. Incubate in ice for 10 min.
- c) After incubation, centrifuge at 12,000 rpm for 2-5 min, carefully transfer the supernatant to a new EP tube and place on ice for later use.
- d) Take a small amount of supernatant (1-2 μL) to measure the protein concentration Cpr using the Bradford method. To ensure the accuracy of the experiment, it is recommended to achieve a protein concentration of 1-4 mg/mL.

3. Caspase enzyme assay:

- 1) Refer to the following table to prepare the groups. This assay can be performed directly in a black 96-well plate. Or transfer the samples to a black 96-well plate before detection.

Sample group (experimental group and untreated control group)	50 μL supernatant (If the volume is less than 50 μL , adjust volume to 50 μL with Cell Lysis Buffer)
Negative control group	50 μL Cell Lysis Buffer

- 2) Add 50 μL of 2X Reaction Buffer (containing DTT) to each well.
- 3) Add 5 μL of 1 mM WEHD-AFC to each well. Incubate at 37°C in the dark for 1-2 h.

***Note:** The negative control group does not need to add VDAD-AFC substrate.

- 4) Measure the Fluorescence signal by microplate reader at Ex/Em = 400/505 nm. Or measure the signal using a fluorometer.

***Note:** This kit detects fluorescence, so a black plate is recommended for microplate readers.

4. Analysis: The following formula can be used to calculate the enzyme activity

$$\text{Caspase activity (\%)} = \frac{F_{\text{Experimental group}} - F_{\text{Negative control group}}}{Cpr_{\text{experimental group}}} \div \frac{F_{\text{untreated control group}} - F_{\text{Negative control group}}}{Cpr_{\text{untreated control group}}} \times 100\%$$

Note: F is the fluorescence signal value for each group; Cpr is the protein concentration; If the sample set is diluted with Cell Lysis Buffer, the F of the original sample needs to be corrected first, that is, $(F_{\text{Experimental group}} - F_{\text{Negative control group}}) \times a$, a is the dilution coefficient of the original sample.

Note

Problems	Cause	Solution
Assay not working	- Cells did not lyse completely - Experiment was not performed at optimal time after apoptosis induction - Plate read at incorrect wavelength - Old DTT used	- Resuspend the cell pellet in the lysis buffer and incubate as described in the datasheet - Perform a time-course induction experiment for apoptosis - Check the wavelength listed in the datasheet and the filter settings of the instrument - Always use freshly thawed DTT in the cell lysis buffer
High Background	- Increased amount of cell lysate used - Increased amounts of components added due to incorrect pipetting - Incubation of cell samples for extended periods - Use of expired kit or improperly stored reagents - Contaminated cells	- Refer to datasheet and use the suggested cell number to prepare lysates - Use calibrated pipettes - Refer to datasheet and incubate for exact times - Always check the expiry date and store the individual components appropriately - Check for bacteria/ yeast/ mycoplasma contamination
Lower signal levels	- Cells did not initiate apoptosis - Very few cells used for analysis - Use of samples stored for a long time - Incorrect setting of the equipment used to read samples - Allowing the reagents to sit for extended times on ice	- Determine the time-point for initiation of apoptosis after induction (time-course experiment) - Refer to datasheet for appropriate cell number - Use fresh samples or aliquot and store and use within one month for the assay - Refer to datasheet and use the recommended filter setting - Always thaw and prepare fresh reaction mix before use
Samples with	Uneven number of cells seeded in the wells	Seed only equal number of healthy cells (correct

erratic readings	• Samples prepared in a different buffer • Adherent cells dislodged and lost at the time of experiment • Cell/ tissue samples were not completely homogenized • Samples used after multiple freeze-thaw cycles • Presence of interfering substance in the sample • Use of old or inappropriately stored samples	- passage number) • Use the cell lysis buffer provided in the kit • Perform experiment gently and in duplicates/ triplicates; apoptotic cells may become floaters • Use Dounce homogenizer (increase the number of strokes); observe efficiency of lysis under microscope • Aliquot and freeze samples, if needed to use multiple times • Troubleshoot as needed • Use fresh samples or store at correct temperatures until use
Unanticipated results	• Measured at incorrect wavelength • Cell samples contain interfering substances	• Check the equipment and the filter setting • Troubleshoot if it interferes with the kit (run proper controls)
General issues	• Improperly thawed components • Incorrect incubation times or temperatures • Incorrect volumes used • Air bubbles formed in the well/tube • Substituting reagents from older kits/ lots • Use of a different 96-well plate	• Thaw all components completely and mix gently before use • Refer to datasheet & verify the correct incubation times and temperatures • Use calibrated pipettes and aliquot correctly • Pipette gently against the wall of the well/tubes • Use fresh components from the same kit • Fluorescence: Black plates; Absorbance: Clear plates

Note: The most probable cause is listed under each section. Causes may overlap with other sections.

APExBIO Technology

www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com