


EdU Cell Proliferation Kit (TMB)

Introduction

Measuring cell proliferation and cell cycle are a fundamental method to assess cell health, determine genotoxicity, and evaluate drug's pharmacodynamic effect. The common method is measuring DNA synthesis directly. In previous experiments, there are several approaches such as the incorporation of radioactive nucleosides (^3H -thymidine) or BrdU. Here, we introduce one new method, click chemistry - CuAAC (Copper-Catalyzed Azide-Alkyne Cycloaddition), and the use of this reaction in direct measurement of S-phase DNA synthesis in cell cycle.

A nucleoside analog of thymidine, EdU (5-ethynyl-2'-deoxyuridine), can be incorporated into DNA strand during DNA synthesis. The alkynyl group of EdU is a biologically inert group that will undergo an extremely selective reaction with dye's azido via a CuAAC reaction to afford an 1,2,3-triazole product. EdU and Biotin azide possess biologically unique moieties to label DNA of proliferating cells, producing low backgrounds and high detection sensitivities. This CuAAC reaction affords superior regioselectivity and quantitative transformation under extremely mild conditions.

EdU Imaging Kits (TMB) specifically labels the DNA of proliferating cells after biotin azide is ligated to EdU, then adds horseradish peroxidase-labeled streptavidin (HRP-Streptavidin) to bind biotin, developing by TMB chromogen, and finally visualized by microplate reader.



Components and Storage

Components	500 Tests	2000 Tests	Storage
EdU 10mM (Component A)	200 µL	800 µL	-20°C
Biotin azide (Component B)	26 µL	110 µL	-20°C away from light
10X EdU Reaction Buffer (Component C)	3.5 mL	13 mL	4°C
CuSO ₄ (100 mM Aqueous Solution) (Component D)	1.2 mL	4.8 mL	4°C
10X EdU Buffer Additive (Component E)	1 vial	4 vials	-20°C
Streptavidin-HRP (Component F)	110 µL	450 µL	-20°C away from light
Streptavidin-HRP Diluent (Component G)	10 mL	40 mL	-20°C
TMB Chromogen Solution (Component H)	50 mL	200 mL	4°C away from light
Shipping: Blue ice		Shelf life: 1 year	

Protocol

1. Preparing reagents

- 1) Allow vials to warm to room temperature before opening.
- 2) Prepare a working solution of 1X EdU reaction buffer (Component C): Dilute Component D in deionized water to make the 1X EdU reaction buffer. After use, store any remaining 1X solution at 2–6°C. When stored as directed, this 1X solution is stable for up to 6 months.
- 3) To make a 10X stock solution of the EdU buffer additive (Component E): Add 0.26 mL of deionized water to the vial and mix until the EdU buffer additive is fully dissolved. After use, store any remaining stock solution at -20°C. When stored as directed, the stock solution is stable for up to 1 year. If the solution develops a brown color, it has degraded and should be discarded.

2. Materials required but not provided

- A. Phosphate-buffered saline (PBS, pH 7.2–7.6)
- B. Fixative (for example, 3.7% Formaldehyde in PBS)
- C. Permeabilization reagent (for example, 0.3% Triton X-100 in PBS)
- D. Endogenous peroxidase blocking solution (for example, 0.3% H₂O₂ in PBS)
- E. Deionized water

3. Labeling cells with EdU (96-well plate)

In initial experiments, we recommend testing a range of EdU concentrations to determine the optimal concentration for your cell type and experimental conditions. If you are currently using a BrdU-based assay for

cell proliferation, a similar concentration to BrdU is a good starting concentration for EdU.

- 1) Plate the cells at the desired density in a 96-well plate, then allow them to recover overnight before additional treatment.
- 2) Prepare a 2X working solution of EdU in a complete medium from the 10 mM stock solution. You can take a final concentration of 10 μ M as a start point. (For example, for a final concentration of 10 μ M, prepare a 2X working solution of 20 μ M).
- 3) Prewarm the 2X EdU solution, then add an equal volume of the 2X EdU solution to the volume of media containing cells to be treated to obtain a 1X EdU solution. (For example, for a final concentration of 10 μ M, replace half of the media with fresh media containing 20 μ M of EdU). We do not recommend replacing all of the media, because this could affect the rate of cell proliferation.
- 4) Incubate the cells for the desired length of time under conditions optimal for the cell type. The time of EdU exposure to the cells allows for direct measurement of cells synthesizing DNA. The choice of time points and the length of time depends on the cell growth rate. Pulse labeling of cells by brief exposures to EdU permits studies of cell-cycle kinetics.
- 5) Incubate under conditions optimal for the cell type for the desired length of time. The time of EdU exposure to the cells allows for the direct measurement of cells synthesizing DNA. The choice of time points and length of time for pulsing depends on the cell growth rate.
- 6) Proceed immediately to Cell fixation and permeabilization.

4. Fixation and permeabilization

- 1) After incubation, remove the media and add 100 μ L of 3.7% formaldehyde in PBS to each well containing the coverslips. Incubate for 15 minutes at room temperature.
- 2) Remove the fixative and wash the cells in each well twice with 100 μ L of PBS.
- 3) Remove the wash solution. Add 100 μ L of 0.3% Triton X-100 in PBS to each well, then incubate at room temperature for 10-15 minutes.
- 4) Wash the cells in each well twice with 100 μ L of PBS.
- 5) Add 100 μ L of 0.3% H₂O₂ PBS per well, and incubate for 10–15 min at room temperature to inactivate endogenous catalase. Subsequently, wash cells 2 times with PBS.

5. Click reaction

- 1) Prepare 1X EdU buffer additive by diluting the 10X stock solution 1:10 in deionized water. Prepare this solution fresh and use the solution on the same day.
- 2) Prepare the Click reaction cocktail. Use this within 15 minutes after preparation.

Components	96-well plate			
	10	20	40	50
1X EdU Reaction Buffer	430 μ L	860 μ L	1.8 mL	2.2 mL
CuSO ₄ (Component D)	20 μ L	40 μ L	80 μ L	100 μ L
Biotin azide (Component B)	0.5 μ L	1 μ L	2 μ L	2.5 μ L
1X EdU Buffer Additive	50 μ L	100 μ L	200 μ L	250 μ L
总体积	500 μ L	1 mL	2 mL	2.5 mL

- 3) Add 50 μ L of click reaction cocktail to each well. (You can adjust the volume of the cocktail to your preference according to previous experiments.) Rock the plate briefly to ensure that the reaction cocktail is distributed evenly.
- 4) Incubate the plate for 30 minutes at room temperature, protected from light (water can be added to the gaps of the plate during incubation to minimize reaction evaporation).
- 5) Remove the reaction cocktail, then wash each well once with 100 μ L of PBS. Remove the wash solution.

6. TMB staining of samples

- 1) In between click reaction incubations, prepare the Streptavidin-HRP working solution regarding the table below and mix thoroughly. Streptavidin-HRP working solution must be prepared and used now and cannot be cryopreserved.

Components	96-well plate			
	10	20	40	50
Streptavidin-HRP (Component F)	2 μ L	4 μ L	8 μ L	10 μ L
Streptavidin-HRP Diluent (Component G)	198 μ L	396 μ L	792 μ L	990 μ L
Total volume	200 μ L	400 μ L	800 μ L	1 mL


- 2) Add 20 μ L Streptavidin-HRP working solution per well and incubate for 30 minutes at room temperature (water can be added to the gaps of the plate during incubation to minimize reaction evaporation).
- 3) After incubation, wash cells three times with PBS, 3-5 min per time.
- 4) Add 100 μ L TMB Chromogen Solution per well, and incubate for 5-30 min at room temperature. During this period, the color can be observed under the microscope several times, and if the color is dark, it can be 5 minutes.
- 5) Absorbance is determined directly at 370 nm or 620-650 nm. Alternatively, add 25 μ L of 2M H₂SO₄ to stop the reaction, and then measure the absorbance at 450 nm.

Note

1. This kit measures absorbance, and it is generally recommended to use a 96-well plate or a 394-well plate for

detection. If you have an instrument that can inspect large plates such as 6-well plates, you can choose according to your habits.

2. If hydroxyurea is required as a control, it can be purchased by our catalog number (B2102 Hydroxyurea).
3. If more EdU is required for animal experiments, it can be purchased by our catalog number (B8337 5-EdU).
4. For your safety and health, please wear lab coats and gloves during the experiment.
5. For research use only. Not to be used in clinical diagnostic or clinical trials.



APExBIO Technology

www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com

