

Cell Membrane Staining Kit (DiO)

Introduction

DiO, also called DiOC₁₈(3), is a green fluorescent probe that is widely used as a lipophilic tracker in living and fixed tissues and cells. DiO can insert into the membrane, and diffuse rapidly, staining the entire membrane surface green. DiO usually does not affect cell viability and development, so it is commonly used in cell tracking, such as anterograde and retrograde neuronal tracking. DiO staining is usually less intense than Dil staining, and occasionally fails completely in fixed tissues.

DiO staining can be used in conjunction with immunofluorescence. It is recommended to use formaldehyde (PFA) for fixation after DiO staining. Meanwhile, if permeabilization is required, it is suggested to use Triton X-100 or digitonin. But permeabilization may affect the localization of DiO in the cell membrane. In addition, DiO has many other applications including detecting cell-cell fusion and adhesion, tracking cell migration, and labeling lipoproteins.

This kit is an optimized DiO staining kit that not only provides ready-to-use DiO dye, but also provides a Staining Enhancer for better staining results and lower staining background. For 96-well plates, the kit can be used 1000 times using 100 µL of staining solution per well.

Composition and storage conditions

| Components | Size | 1000 Assays | Storage |
|--------------------------|------|--------------------|-----------------------|
| DiO (200×) | | 0.5 mL | -20°C away from light |
| Staining Enhancer (200×) | | 0.5 mL | -20°C |
| Staining Buffer | | 100 mL | -20°C |
| Shipping: Blue ice | | Shelf life: 1 year | |

Experimental manipulation

1. Preparation of the working solution

- 1) Refer to the following table to prepare the working solution, the working solution is unstable and needs to be prepared freshly. DiO (200×) and Staining Enhancer (200×) should be thoroughly mixed before use. For 96-well plates, 100 µL of staining solution is required per well. For 6-well plates, 1 mL of staining solution is required per well.

| | | | |
|--------------------------|--------|---------|--------|
| Number of samples | 10 | 100 | 1000 |
| DiO (200×) | 5 µL | 50 µL | 500 µL |
| Staining Enhancer (200×) | 5 µL | 50 µL | 500 µL |
| Staining Buffer | 990 µL | 9.90 mL | 99 mL |
| Total | 1 mL | 10 mL | 100 mL |

***Note:** The working solution concentration can be adjusted and optimized according to the cells and experimental system.

2. Staining of suspension cells

- 1) Harvest cells and centrifuge at 1000 rpm for 5 min, remove the supernatant. Suspend the cell pellet in a suitable working solution at a density of $1-2 \times 10^6$ cells/mL.
- 2) Incubate at 37°C away from light for 2-20 min. The optimal incubation time varies depending on the cell type.
- 3) Centrifuge at 1000 rpm for 5 min, remove the supernatant.
- 4) Gently resuspend the cell pellet in a pre-warmed growth medium, and wash 2 times.
- 5) Gently resuspend the cell pellet in a suitable buffer (serum-free medium or PBS), then watch by a microscope or flow cytometry (Ex/Em= 484/501 nm). DiO can be watched directly with a FITC filter.

3. Staining of adherent cells

- 1) Grow adherent cells on cell dishes, cell plates or sterile glass coverslip.
- 2) Remove the growth medium and wash cells twice with PBS.
- 3) Add the working solution to cover the cells. For 96-well plates, 100 µL of staining solution is required per well. For 6-well plates, 1 mL of staining solution is required per well. Incubate at 37°C away from light for 2-20 min. The optimal incubation time varies depending on the cell type.
- 4) Remove the working solution and wash the cells twice with PBS.
- 5) Gently add a suitable buffer (serum-free medium or PBS) to cover the cells, then watch by a microscope (Ex/Em= 484/501 nm). DiO can be watched directly with a FITC filter.

4. Fix after staining

- 1) If further immunofluorescence experiments are required after staining, fix cells with 4% PFA.
- 2) If permeabilization is also required, permeabilization with 0.1% Triton X-100 or digitonin is recommended. However, permeabilization is likely to affect the staining of cell membranes by DiO and increase the staining background.
- 3) In addition, detergents may dissolve lipids on the cell membrane and affect the localization of the cell membrane of DiO, so blocking solutions, antibody diluents, and washes should not contain detergents. It is recommended to mount directly with PBS. Do not use mounting media containing glycerin or other

organic matter, as this will affect the staining effect.

5. Staining after fixation

- 1) Fix the cells with 4% PFA.
- 2) If permeabilization is also required, permeabilization with 0.1% Triton X-100 or digitonin is recommended. However, permeabilization is likely to affect the staining of cell membranes by DiO and increase the staining background.
- 3) Perform the immunofluorescence staining as usual. Detergents may dissolve lipids on the cell membrane and affect the membrane localization of DiO, so blocking solutions, antibody diluents, and washes should not contain detergents.
- 4) Add the appropriate working solution to cover the cells. For 96-well plates, 100 μ L of staining solution is required per well. For 6-well plates, 1 mL of staining solution is required per well. Incubate at 37°C away from light for 2-20 min. The optimal incubation time varies depending on the cell type.
- 5) Remove the working solution and wash the cells twice with PBS.
- 6) Gently add a suitable buffer (serum-free medium or PBS) to cover the cells, then watch by a microscope (Ex/Em= 484/501 nm). DiO can be watched directly with a FITC filter.

Note

1. Fluorescent probes are easy to quench, please protect them from light when using.
2. If the labeling time is too long or cells are cultured after staining, the probe may also enter the cell to stain other organelles.
3. For your safety and health, please wear lab coats and gloves during the experiment.
4. For research use only. Not to be used in clinical diagnostic or clinical trials.



APEx BIO Technology
www.apexbt.com
7505 Fannin street, Suite 410, Houston, TX 77054.
Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com