

ELISA Coating Buffer (10×)

Introduction


This product is a concentrated solution with a pH of 9.6. It must be diluted with deionized water to a 1× concentration before use. It is used for coating antigens or antibodies onto ELISA microplates in enzyme-linked immunosorbent assay (ELISA) procedures. In practice, simply take an appropriate amount of this product as the coating buffer and dilute the antigen or antibody solution to the desired coating concentration. Then, add 100 µL of the diluted solution to each well of the ELISA microplate, and incubate overnight at 4°C for coating.

Protocol

1. Dilute to 1× with deionized water before use.
2. Take an appropriate amount of the coating buffer to dilute the antigen or antibody solution to the corresponding coating concentration. Add 100 µL to each well of the ELISA microplate and incubate overnight at 4°C.

Note

1. Storage condition: Store at room temperature; stable for one year.
2. This product is for scientific use only.



APExBIO Technology
www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com

