

# HyperChrom Protein G 4FF Agarose

## Product description

HyperChrom Protein G 4FF Agarose is an affinity chromatography medium made by coupling Protein G to the Chromstar substrate by CNBr activation. Protein G can bind to the Fc region of a variety of mammalian IgG and has a weak affinity for its Fab region.

Compared to Protein A, Protein G has a broader binding spectrum. Therefore, HyperChrom Protein G 4FF Agarose can be used to isolate and purify immunoglobulins and their fragments from different liquid samples such as cell culture media or serum. HyperChrom Protein G 4FF Agarose can also be used to isolate immune complexes.

## Components and storage conditions

100 AS	Coore.	
Components	PC2006-25 mL	PC2006-100 mL
HyperChrom Protein G 4FF Agarose	25 mL	100 mL
Store the components at 1°C for 5 years		

Store the components at 4°C for 5 years.

# Product parameters

HyperChrom Protein G 4FF Agarose chromatography media parameters

Name	Description	
Chromatography media type	Affinity chromatography media	
Ligation	Protein G (E. coli recombinant expression).	
Scaffolding	Highly cross-linked agarose	
Average particle size	90 μm	
Ligand density	~2 mg Protein G/mL chromatography medium	
Dynamic load	$\geq$ 20 mg human IgG/mL chromatography medium*	
Flow rates are recommended	90-300 cm/h	
Maximum flow rate	400 cm/h	
Withstand pressure	0.3 MPa	
Use temperature	4-30°C	
pH stability **	3-10	
Chemical resistance	Common aqueous solution, 30% isopropanol*** 75% ethanol*** 6 M	
	guanidine hydrochloride, 8 M urea	
*Dynamic load measurement conditions: column loading height: 10 cm, test flow rate 200 cm/h, test buffer:		
0.02 M NaH <sub>2</sub> PO <sub>4</sub> solution, pH 7.0, test sample: 2 mg/m IgG sample, when the penetration of IgG reaches		

APENBIC

10%, unit medium volume (mL) IgG loading volume (mg).

\*\*After 7 days of storage in the environment of 40 °C and pH 3-9, the physicochemical properties and functions of the chromatography medium did not change significantly.

\*\*\* 30% and 75% are the volume ratio (v/v).

## **Experimental manipulation**

Experimental manipulation 1. Preparation of buffers	APPER Concernent
Buffer type	Buffer components
Balance/Bind/Wash Buffer	0.02-0.05 M PB or Tris, pH 7.0-8.0, 0.15 M NaCl can be added to reduce non-specific adsorption.
Elution Buffer	0.1 M glycine, pH 2.5-3.0

**Sample preparation** 2.

Prepare samples for purification.

#### **Chromatographic conditions** 3.



- Flow rate selection: Linear flow rate of 90-300 cm/h is generally selected according to the height of the column bed.
- Sample preparation: To prevent the sample from clogging the column, the sample needs to be filtered with a  $0.2/0.45 \,\mu\text{m}$  microporous membrane (after inclusion body disruption) before loading, and it is recommended that the pH and conductivity of the sample be adjusted to be consistent with the equilibration buffer (the pH and conductivity of the sample can be adjusted by dilution, ultrafiltration, and desalting).

#### Loading columns 4.

The following column loading methods are suitable for filling laboratory-scale chromatography columns:

- 4.1 Supplies required for column mounting.
- (1) Chromatographic medium: HyperChrom Protein G 4FF Agarose.
- (2) Chromatographic empty column: laboratory-scale chromatography empty column and column loader.
- (3) Solution required:
  - Column loading solution: 20% ethanol. a)
  - Exhaust solution: 20% ethanol. b)
- (4) Column loading tools: sand core funnel, stirring rod, measuring cylinder, etc.

- 4.2 Preparation before column loading
- Calculate the volume of chromatography medium Vm required for column loading (volume of the chromatography medium part after sufficient sedimentation), and calculate the formula:

Vm = cross-sectional area of the chromatography column X height of the column bed where the column is planned to be loaded x compression ratio of the chromatography medium.

(Note: HyperChrom Protein G 4FF Agarose has a compression ratio of 1.15).

- (2) The chromatography medium was transferred to the sand core funnel, and the column loading solution was cleaned and filtered with about 3 times the volume of the chromatography medium, and the column chromatography medium to be loaded was replaced with the column loading solution.
- (3) For the preparation of the glue suspension of the column chromatography medium to be loaded, the suitable proportion of column loading glue suspension for HyperChrom Protein G 4FF Agarose chromatography medium is 45%-55%. To obtain an accurate chromatography medium volume, the chromatography medium can be placed in a graduated cylinder and settled overnight or centrifuged at low speed (3000 rpm, 5 min) to simulate the natural sedimentation effect of the chromatography medium, and then measured.
- (4) Check the empty column to be used to ensure it is clean and leak-free.
- 4.3 Column mounting
- (1) Vent the column bottom membrane (screen) with 20% ethanol.
- (2) After sufficient exhaust, screw the plug or close the column bottom valve at the bottom interface of the column, and continue to inject a small amount of 20% ethanol until the bottom of the column is covered.
- (3) Adjust the chromatography column to vertical.
- (4) The column head is connected to the chromatography system, through which a low flow rate of 5 m L/min is provided, and the column head filter membrane (screen) is vented with 20% ethanol.
- (5) Thoroughly stir the prepared chromatography medium suspension with a stir bar, and then slowly pour into the prepared chromatography empty column at one time.

Note: If the volume of the glue suspension exceeds the empty column volume, it should be extended by using a column loader or connecting another empty column tube with a connector.

- (6) Place the gaseated stigma into the chromatography column, fully fit the glue suspension level, and remove all air bubbles. Then tighten the column head seal.
- (7) Start the system pump, adjust the flow rate to 300 cm/h, and use the liquid flow to press the column bed. During this period, the pressure should not exceed 0.3 MPa. If the pressure is over, the flow rate

CID 26 mm 50 mm 10 mm 16 mm V/FR LFR 5.3 mL/ min 19.6 mL/min 60 cm/h 0.8 mL/min 2.0 mL/min 8.8 mL/min 100 cm/h 1.3 mL/min 3.3 mL/min 32.7 mL/min 1 50 cm/h 2.0 mL/min5.0 mL/min 13.3 mL/ min 49.1 mL/min 200 cm/h 2.6 mL/min6.7 mL/min 17.7 mL/ min 65.4 mL/min 3.9 mL/min 300 cm/h 10.0 mL/min | 26.5 mL/min 98.1 mL/min 600 cm/h 7.9 mL/min 20.1 mL/ min 53.1 mL/min 196.3 mL/ min Note:

needs to be reduced (see the table below for flow rate conversion).

CID: Chromatographic inner diameter

V/FR: Volumetric flow rate

LFR: Linear flow rate

Table 1 Flow rate conversion table of different specifications of chromatography columns

- (8) After the column bed is stabilized (the glue surface no longer falls), mark the position of the glue surface at this time. Stop the pump and press the column head down to 3-5 mm below the marked position.
- (9) Re-apply the flow rate of 600 cm/h, if the glue surface no longer drops, that is, the column loading is completed. If the glue surface falls, repeat steps 8-9.

Note: The recommended workflow speed does not exceed 75% of the flow rate of the loading column.

#### Column efficiency determination (optional). 5.

Select one of the two test methods shown in the table below for column effectiveness testing. Use the mobile phase equilibrium chromatography column to the baseline to be stable, load the sample into the chromatography column, continue to use the mobile phase for rinsing, and after the chromatographic peak is completed to return to the baseline, end the run, integrate the chromatographic peak, and evaluate the loading effect.

Table 2 Statistical table of two column efficiency measurement methods

	Acetone method	NaCl method
sample	1% (v/v) acetone in water	2 M NaCl in water
Sample volume	1% column volume	1% column volume
Mobile phase	water	0.2 M NaCl in water
velocity of flow	30 cm/h	30 cm/h
Detector	UV 280 nm	electrical conductivity

The main evaluation criteria for the effect of column loading are N/m (number of plates per meter) and As (symmetry factor), which are calculated as follows:

$$\frac{N}{m} = 5.54x (\frac{V_R}{W_h})^2 \times \frac{1}{L}$$

$$As = b/a$$

Column efficiency qualification standards: N/m > 3000, 0.8 < As < 1.5

\*Parameter Notes:

L = column height, VR = reserved volume, Wh = half-peak width, a = left half-peak width at 10% peak height, b = right half-peak width at 10% peak height

#### 6. Chromatographic steps

- Equilibrium: Use Balance/Bind/Wash Buffer to fully equilibrate the column to pH and conductivity stable and substantially consistent with the equilibration buffer, which typically requires 3-5 times the column volume.
- (2) Sample loading: Determine the sample loading volume and loading volume on the HyperChrom Protein G 4FF Agarose based on the binding load measured in the pilot experiment.
- (3) Washing: Rinse the column with Balance/Bind/Wash Buffer or other suitable buffers until UV stable and return to baseline.
- (4) Elution: Elution Buffer is lowered for washing. Immediately after elution, the collection components should be neutralized with 1 M Tris-HCl, pH 9 buffer, which can be added at a ratio of 60-200 µl to 1 ml of the elution collection.
- (5) Re-equilibration: Re-equilibrate the chromatography column with a 10x volume Balance/Bind/Wash Buffer.

### 7. Cleaning and recycling



As the number of uses of the chromatography medium increases, contaminants (e.g., lipids, endotoxins, proteins, etc.) accumulate on the chromatography column. Regular in-place cleaning is essential to keep the column in stable working condition. Determine the frequency of in-place cleaning according to the degree of contamination of the chromatography medium (if the contamination is serious, it is recommended that in-place cleaning should be carried out after each use to ensure repeatable results and extend the working life of the chromatography medium).

For different types of impurities and contaminants, cleaning can be carried out under the following conditions:

- Method 1: To remove strong hydrophobic proteins, lipoproteins and lipids, etc., it can be treated with 0.1% nonionic detergent at 37 °C with a contact time of 1 min, and then rinsed with at least 5 column volumes of conjugate solution.
- Method 2: Soak in 70% ethanol for 12 h to remove lipids, then rinse with at least 5 column volumes of binding solution.

### 8. sterilization

To reduce the microbial load, 20% ethanol was used for more than 6 h to achieve sterilization.

#### 9. stockpile

Unopened chromatography media, please keep in the original container; The completed chromatography column is first soaked with 20% ethanol solution and then the upper and lower column heads are closed. The storage environment is  $2 \sim 8$  °C.

#### **10.** Destruction and recycling

- Since HyperChrom Protein G 4FF Agarose chromatography media is difficult to degrade in nature, incineration of discarded chromatography media is recommended to protect the environment.
- For chromatography media exposed to bioactive samples such as viruses and blood, please follow local biosafety requirements before destroying or disposing of them.

### Notes

- It is recommended that the buffer and protein solution used for purification be filtered through a 0.22 μm or 0.45 μm membrane and then used on the column.
- 2. This product is for scientific purposes only.



APENER BIO









### APExBIO Technology www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054. Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com











