

Cat. No. K1094

T5 Exonuclease (RNase-free)

Introduce

T5 Exonuclease is a nucleic acid exonuclease that degrades double- or single-stranded DNA in a 5' \rightarrow 3' direction. Degradation can begin either from the 5' end of double- or single-stranded DNA, or from gaps or nicks in linear or cyclic double-stranded DNA, but not from superhelix double-stranded DNA. T5 Exonuclease is often used for Gibson Assembly; Degradation of linear single-stranded, double-stranded DNA or absent plasmid DNA; Removal of incomplete ligation products from the linked cyclic double-stranded DNA; Degradation of linear and absent plasmid DNA to obtain high-purity supercoiled plasmid DNA; Removal of metamorphic DNA produced during plasmid extraction by alkali cleavage; Improve transfection efficiency of small extract plasmid cDNA libraries.

Composition and storage conditions

| Components Size | 1000 U | 5000 U | 10000 U |
|------------------------------------|---------------|---------|---------|
| T5 Exonuclease (RNase-free) | 0.1 ml | 0.5 ml | 1 ml |
| 10X T5 Exonuclease Reaction Buffer | 1 ml | 1 ml x5 | 1ml x10 |

Experimental operation

 1. Configure the reaction system in the ice bath according to the following table:

 Reagent
 Volume

Final Contract of the following table:

| Reagent | Volume | Final Concentration |
|---------------------|---------|---------------------|
| DNA | Xμl | 1 µg |
| 10X Reaction buffer | 5µl | 1X |
| T5 Exonuclease | 1µl | 10 U |
| Nuclease-free H2O | to 50µ1 | - |

- 2. Gently mix the reaction system up and down with a pipette and centrifuge quickly.
- 3. Incubate at 37 °C for 10~30min.
- 4. Stop the reaction with 6X Purple Gel Loading Dye or adding EDTA to at least 11 mM.
- 5. The digested products are analyzed by agarose gel or polypropylene gel electrophoresis.

Note: T5 Exonuclease should be added to the reaction system last, and the mixing reaction system should be paid attention

to before joining; Store in an ice box or on an ice bath.

Note

- 1. Definition of Enzyme Activity Units: The amount of enzyme required to catalyze the production of 1 nmol of DNA from a double-stranded DNA substrate within 30 minutes in a 50 μ L reaction system at 37 °C.
- Liquid components in T5 nucleic acid-free exonuclease: 50 mM Tris-HCl (pH 7.5, 25 °C), 100 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.1% TritonX-100, 50% glycerol.
- 3. Buffer components in the 10X T5 Exonuclease Reaction Buffer: 500 mM Potassium Acetate, 200 mM Tris-acetate (pH 7.9, 25°C), 100 mM Magnesium Acetate, 10 mM DTT.
- 4. T5 Exonuclease (RNase-free) quality and experimental validation:
 - I. Protease purity > 95%.
 - II. Verification of activity:
 - 10U of T5 Exonuclease and 200ng of 300 base RNA substrate reacted at 37 °C for 4 h with no change in the electrophoresis band of RNA.
 - 50U of T5 Exonuclease and 1 µg of supercoiled PhiX174 DNA reacted at 37 °C for 4 h without alteration of the electrophoretic band of the DNA.

5. T5 Exonuclease inactivation conditions: T5 Exonuclease can be inactivated by adding a final concentration greater than 11 mM EDTA to terminate the reaction, or a DNA loading buffer containing SDS (0.08% final concentration of SDS).

APENBIO MARCA

APENEN BIO



