



APEX BIO
Achieve Perfection, Explore the Unknown

High Performance Liquid Chromatography


1. Analysis Information

Product Name:	Scriptaid
Operator:	David
Injection Date:	11/12/2013 7:09:50 PM
Batch No.:	1

2. HPLC Condition

Column:	Waters XBridge C18, 3 μ m, 3 mm \times 100 mm																
Solvent A:	0.1% H3PO4 in 100% Acetonitrile																
Solvent B:	0.1% H3PO4 in 100% Water																
Gradient:	<table border="1"> <thead> <tr> <th>Time</th> <th>A</th> <th>B</th> </tr> </thead> <tbody> <tr> <td>0.0min</td> <td>95%</td> <td>5%</td> </tr> <tr> <td>5.0min</td> <td>95%</td> <td>5%</td> </tr> <tr> <td>10.0min</td> <td>95%</td> <td>5%</td> </tr> <tr> <td>10.0min</td> <td>Stop</td> <td>Stop</td> </tr> </tbody> </table>		Time	A	B	0.0min	95%	5%	5.0min	95%	5%	10.0min	95%	5%	10.0min	Stop	Stop
Time	A	B															
0.0min	95%	5%															
5.0min	95%	5%															
10.0min	95%	5%															
10.0min	Stop	Stop															
Flow rate:	0.5ml/min																
Wavelength:	254nm																
Volume:	10 μ l																

3. Result



Rank	RetTime(min)	Area(mAU*s)	Area %
1	3.803	23.05198	0.3026
2	5.142	7579.52588	99.4910
3	6.453	15.72489	0.2064