

Product Information

Caspase-9 Fluorometric Assay Kit

I. Kit Contents:

Component	K2018-25 25 assays	K2018-100 100 assays	K2018-200 200 assays	K2018-400 400 assays	Part Number
Cell Lysis Buffer	25 ml	100 ml	100 ml	100 ml	K2018-C-1
2X Reaction Buffer	2 ml	4 x 2 ml	16 ml	32 ml	K2018-C-2
LEHD-AFC (1 mM)	125 μ l	500 μ l	2 x 0.5 ml	2 x 1 ml	K2018-C-3
DTT (1 M)	100 μ l	400 μ l	400 μ l	400 μ l	K2018-C-4

II. Introduction:

Cysteine-dependent aspartate-directed proteases (Caspases) are a family of cysteine proteases that play important roles in apoptosis, necrosis, and inflammation. Sequential activation of caspases plays an important role in cell apoptosis. Caspase-9 is a member of the caspase family. Caspase-9, as well as caspase-2, -8 and -10, is an initiator caspase. Caspase-9 is involved in the mitochondrial death pathway and is activated during apoptosis. The activation of JNK/SAPK stress signaling pathways induces the release of cytochrome c from mitochondria and activation of apaf-1, which then cleaves caspase-9 into the active form. The active Caspase-9 cleaves and activates Caspase-3 and Caspase-7, which then cleave poly ADP ribose polymerase.

Caspase-9 Fluorometric Assay Kit provides a convenient and simple way for detecting the LEHD-dependent caspase activity. LEHD-AFC (AFC: 7-amino-4-trifluoromethyl coumarin) emits blue light ($\lambda_{max} = 400$ nm); while cleavage of the substrate LEHD-AFC by Caspase-9 or related caspases, free AFC emits a yellow-green fluorescence ($\lambda_{max} = 505$ nm), which can be quantified by using a fluorescence microtiter plate reader or a fluorometer. Comparison of the fluorescence of free AFC from an apoptotic sample with an uninduced control determines the fold increase in Caspase-9 activity.

III. Caspase-9 Assay Protocol:

A. General Considerations

Aliquot enough 2X Reaction Buffer for the number of assays to be performed. Add DTT to the 2X Reaction Buffer immediately before use (10 mM final concentration: add 10 μ l of 1.0 M DTT stock per 1 ml of 2X Reaction Buffer).

Protect LEHD-AFC from light.

B. Assay Procedure

1. Induce apoptosis in cells by desired method. Concurrently incubate a control culture without induction.
2. Count cells and pellet $1 - 5 \times 10^6$ cells or use 50-200 μ g cell lysates if protein concentration has been measured.
3. Resuspend cells in 50 μ l of chilled Cell Lysis Buffer. Incubate cells on ice for 10 minutes.
4. Add 50 μ l of 2X Reaction Buffer (containing 10 mM DTT) to each sample. Add 5 μ l of the 1 mM LEHD-AFC substrate (50 μ M final concentration) and incubate at 37°C for 1-2 hour.
5. Read samples in a fluorometer equipped with a 400 nm excitation filter and 505 nm emission filter. For a plate-reading set-up, transfer the samples to a 96-well plate.

You may also perform the entire assay directly in a 96-well plate.

Fold-increase in Caspase-9 activity can be determined by comparing these results with the level of the uninduced control.

IV. Storage and Stability:

Store kit at -20°C (Store Cell Lysis Buffer and 2X Reaction Buffer at 4°C after opening). All reagents are stable for 1 year under proper storage conditions.

General Troubleshooting Guide for Caspase Colorimetric and Fluorometric Kits:

Problems	Cause	Solution
Assay not working	<ul style="list-style-type: none"> • Cells did not lyse completely • Experiment was not performed at optimal time after apoptosis induction • Plate read at incorrect wavelength • Old DTT used 	<ul style="list-style-type: none"> • Resuspend the cell pellet in the lysis buffer and incubate as described in the datasheet • Perform a time-course induction experiment for apoptosis • Check the wavelength listed in the datasheet and the filter settings of the instrument • Always use freshly thawed DTT in the cell lysis buffer
High Background	<ul style="list-style-type: none"> • Increased amount of cell lysate used • Increased amounts of components added due to incorrect pipetting • Incubation of cell samples for extended periods • Use of expired kit or improperly stored reagents • Contaminated cells 	<ul style="list-style-type: none"> • Refer to datasheet and use the suggested cell number to prepare lysates • Use calibrated pipettes • Refer to datasheet and incubate for exact times • Always check the expiry date and store the individual components appropriately • Check for bacteria/ yeast/ mycoplasma contamination
Lower signal levels	<ul style="list-style-type: none"> • Cells did not initiate apoptosis • Very few cells used for analysis • Use of samples stored for a long time • Incorrect setting of the equipment used to read samples • Allowing the reagents to sit for extended times on ice 	<ul style="list-style-type: none"> • Determine the time-point for initiation of apoptosis after induction (time-course experiment) • Refer to datasheet for appropriate cell number • Use fresh samples or aliquot and store and use within one month for the assay • Refer to datasheet and use the recommended filter setting • Always thaw and prepare fresh reaction mix before use
Samples with erratic readings	<ul style="list-style-type: none"> • Uneven number of cells seeded in the wells • Samples prepared in a different buffer • Adherent cells dislodged and lost at the time of experiment • Cell/ tissue samples were not completely homogenized • Samples used after multiple freeze-thaw cycles • Presence of interfering substance in the sample • Use of old or inappropriately stored samples 	<ul style="list-style-type: none"> • Seed only equal number of healthy cells (correct passage number) • Use the cell lysis buffer provided in the kit • Perform experiment gently and in duplicates/triplicates; apoptotic cells may become floaters • Use Dounce homogenizer (increase the number of strokes); observe efficiency of lysis under microscope • Aliquot and freeze samples, if needed to use multiple times • Troubleshoot as needed • Use fresh samples or store at correct temperatures until use
Unanticipated results	<ul style="list-style-type: none"> • Measured at incorrect wavelength • Cell samples contain interfering substances 	<ul style="list-style-type: none"> • Check the equipment and the filter setting • Troubleshoot if it interferes with the kit (run proper controls)
General issues	<ul style="list-style-type: none"> • Improperly thawed components 	<ul style="list-style-type: none"> • Thaw all components completely and mix gently before use

	<ul style="list-style-type: none"> • Incorrect incubation times or temperatures • Incorrect volumes used • Air bubbles formed in the well/tube • Substituting reagents from older kits/ lots • Use of a different 96-well plate 	<ul style="list-style-type: none"> • Refer to datasheet & verify the correct incubation times and temperatures • Use calibrated pipettes and aliquot correctly • Pipette gently against the wall of the well/tubes • Use fresh components from the same kit • Fluorescence: Black plates; Absorbance: Clear plates
<p>Note# The most probable cause is listed under each section. Causes may overlap with other sections.</p>		

For research use only! Not to be used in humans.

Our promise

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For more details, please visit <http://www.apexbt.com/> or contact our technical team.

Tel: +1-(832)696-8203

Fax: +1-832-641-3177

Email: sales@apexbt.com