

Brilliant Green Stain Solution (1%)

Introduction

Brilliant Green, also known as Basic Green or Brilliant Green SF, is a basic dye. The staining intensity and results of the 1% Brilliant Green (BG) Stain vary depending on the tissue type and application. It is commonly used for staining bacterial specimens.

Storage


Store at room temperature protected from light, stable for 1 year.

Protocol

1. Follow the specific requirements of your experiment.
2. Typical staining time is 1-3 min, which may vary depending on the type of tissue being stained.

Note

1. For your safety and health, please wear lab coats and gloves during the experiment.
2. For research use only. Not to be used in clinical diagnostic or clinical trials.



APExBIO Technology
www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com

