

Tissue Dissociation Solution

Introduction

Tissue Dissociation Solution can efficiently dissociate tissue samples into cell clumps or cell suspensions while maintaining cell viability, which can be used for organoid culture.

Components and Storage

Components	Size	100 mL	Storage
Tissue Dissociation Solution		100 mL	-20°C
Shipping: Dry ice		Shelf life: 6 months	

Protocol

1. The dissociation step varies from tissue to tissue, please refer to the steps in our corresponding organoid culture kit or organoid growth medium protocol.

Note

1. For your safety and health, please wear lab coats and gloves during the experiment.
2. For research use only. Not to be used in clinical diagnostic or clinical trials.

APExBIO Technology
www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com