

Organoid Dissociation Solution

Introduction

Organoid Dissociation Solution is effective for recovering organoids from Matrigel and digesting organoids into small cell clusters while maintaining cell viability, making it ideal for organoid passaging and expansion.

Components and Storage

Components Size	100 mL	Storage
Organoid Dissociation Solution	100 mL	-20°C
Shipping: Dry ice	Shelf life: 6 months	Contraction of the second seco

Protocol

1. The dissociation step varies from organoids, please refer to the steps in our corresponding organoid culture kit or organoid growth medium protocol.

Note

- 1. For your safety and health, please wear lab coats and gloves during the experiment.
- 2. For research use only. Not to be used in clinical diagnostic or clinical trials.

