

## Anti-Adherence Rinsing Solution

### Introduction

Anti-Adherence Rinsing Solution is designed for pre-treating cultureware to effectively reduce the non-specific adhesion of cells and organoids during the culture.

### Components and Storage

Components	Size	100 mL	Storage
Anti-Adherence Rinsing Solution		100 mL	4°C
Shipping: Blue ice		Shelf life: 1 year	

### Protocol

1. Rinse cultureware such as centrifuge tubes, cell culture plates, pipettes, and pipette tips with Anti-Adherence Rinsing Solution in advance before use.

### Note

1. For your safety and health, please wear lab coats and gloves during the experiment.
2. For research use only. Not to be used in clinical diagnostic or clinical trials.

**APExBIO Technology**

**[www.apexbt.com](http://www.apexbt.com)**

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: [info@apexbt.com](mailto:info@apexbt.com)