

Recombinant Human Pleiotrophin

Information

Gene ID	5764		
Accession #	P21246		
Alternate Names			
Source	Escherichia coli.		
M.Wt	Approximately 15.3 kDa, a single non-glycosylated polypeptide chain containin 136 amino acids.		
AA Sequence	GKKEKPEKKV KKSDCGEWQW SVCVPTSGDC GLGTREGTRT GAECKQTMKT QRCKIPCNWK KQFGAECKYQ FQAWGECDLN TALKTRTGSL KRALHNAECQ KTVTISKPCG KLTKPKPQAE SKKKKKEGKK QEKMLD		
Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.		
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles - 12 months from date of receipt, -20 to -70 °C as supplied - 1 month, 2 to 8 °C under sterile conditions after reconstitution - 3 months, -20 to -70 °C under sterile conditions after reconstitution		
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS, pH 7.4.		
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at \leq -20 °C. Further dilutions should be made in appropriate buffered solutions.		
Biological Activity	Fully biologically active when compared to standard. The biological activity was measured by its ability to enhance neurite outgrowth of E16-E18 rat embryonic cortical neurons, when neurons were plated on 96 well culture plates that had been pre-coated with 100 μl/well of a solution of 5-10 μg/ml rHuPTN.		
Shipping Condition	Gel pack.		
Handling	Centrifuge the vial prior to opening.		
Usage	For Research Use Only! Not to be used in humans.		

Components and Storage

Components	5µg	100µg	500µg
Recombinant Human Pleiotrophin	5µg	100µg	500µg

Use a manual defrost freezer and avoid repeated freeze-thaw cycles

- 12 months from date of receipt, -20 to -70 °C as supplied
- 1 month, 2 to 8 °C under sterile conditions after reconstitution
- 3 months, -20 to -70 °C under sterile conditions after reconstitution

Quality Control

Quality Control		C. Control	
Purity	> 96 % by SDS-PAGE and HPLC analyses.	P. Laure tre	
Endotoxin	Less than 0.1 EU/μg of rHuPTN as determine	an 0.1 EU/μg of rHuPTN as determin <mark>ed b</mark> y LAL method.	

Description

Reference









APExBIO Technology

www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com