



Information

Gene ID	246250	
Accession #	Q8K4J7	
Alternate Names		
Source	Escherichia coli.	
M.Wt	Approximately 20.2 kDa, a disulfide-linked homodimer consisting of two 95 amino acid polypeptide chains.	
AA Sequence	MPSMSLCPMD EAISKKINQD FSSLLPAAMK NTVLHCWSVS SRGRLASCPE GTTVTSCSCG SGCGSWDVRE DTMCHCQCGS IDWTAARCCT LRVGS	
Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.	
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles - 12 months from date of receipt, -20 to -70 °C as supplied - 1 month, 2 to 8 °C under sterile conditions after reconstitution - 3 months, -20 to -70 °C under sterile conditions after reconstitution	
Formulation	Lyophilized from a 0.2 μ m filtered concentrated solution in PBS, pH7.4, with 0.02 % Tween-20.	
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at ≤ -20 °C. Further dilutions should be made in appropriate buffered solutions.	
Biological Activity	Test in process.	
Shipping Condition	Gel pack.	
Handling	Centrifuge the vial prior to opening.	
Usage	For Research Use Only! Not to be used in humans.	

Components and Storage

Components	5µg	100µg	500µg
Recombinant Rat Resistin	5µg	100µg	500µg

Use a manual defrost freezer and avoid repeated freeze-thaw cycles

- 12 months from date of receipt, -20 to -70 °C as supplied
- 1 month, 2 to 8 °C under sterile conditions after reconstitution
- 3 months, -20 to -70 °C under sterile conditions after reconstitution

Quality Control

	A STATE OF THE STA	akno.
Purity	> 97 % by SDS-PAGE and HPLC analyses.	P and the state of
Endotoxin	Less than 0.1 EU/µg of rRtResistin as determ	ined by LAL method.

Description

Reference





10





APExBIO Technology

www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com