

Mass Spectrometry


Analysis Information

Product Name:	Cy3 NHS ester (non-sulfonated)
Operator:	David
Injection Date:	3/12/2019 1:15:40 PM
Inj. Volume:	1
Batch No.:	2

Acquisition Parameter

Probe:	ESI
Gas Temp:	350°C
Gas Flow:	11L/min
Nebulizer:	50psi
VCap:	4000v
Fragmentor:	70v
Reference:	m/z 121.050873, 922.009798
T.Flow:	0.3ml/min
B.conc:	20%H2O/80%ACN

Result



Molecular Weight: 590.15

m/z([M-Cl]⁻) Base Peak: 554.3154

Caution

FOR RESEARCH PURPOSES ONLY.

NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

APEX BIO Technology

www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com