

High Performance Liquid Chromatography

Analysis Information

Product Name: Pam3CSK4 TFA

Operator: David

Injection Date: 12/19/2024 16:43:57

Batch No.: 2

HPLC Condition

Column : Athena C18, 3 μ m, 2.1mm $\times$ 100 mm

Solvent A: H2O+0.1%TFA

Solvent B: ACN+0.1%TFA

Solvent C: ME OH

Gradient:


Time	A	B	C
0.0min	30%	0%	70%
2.0min	30%	0%	70%
12.0min	0%	0%	100%
14.0min	0%	0%	100%
18.0min	Stop	Stop	Stop

Flow rate: 1ml/min

Wavelength: 210nm

Volume: 50 μ l

Result



Rank	RetTime(min)	Area(mAU*s)	Area %
1	10.811	18.55	0.60
2	10.943	53.32	1.71
3	12.907	3041.87	97.69

Caution

FOR RESEARCH PURPOSES ONLY.

NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

APExBIO Technology

www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com