

## PBS (Phosphate-Buffered Saline)

### Introduction


PBS (Phosphate-Buffered Saline), containing 135 mM NaCl, 4.7 mM KCl, 10 mM Na<sub>2</sub>HPO<sub>4</sub>, and 2 mM NaH<sub>2</sub>PO<sub>4</sub>, with a pH of 7.3±0.1. This sterile solution is ready for direct use in cell washing or other routine laboratory applications.

### Storage

Store at 4°C for 6 months, or -20°C for 1 year.

### Note

1. Avoid microbial contamination during use.
2. For your safety and health, please wear lab coats and gloves during the experiment.
3. For research use only. Not to be used in clinical diagnostic or clinical trials.



**APEX BIO Technology**  
[www.apexbt.com](http://www.apexbt.com)

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: [info@apexbt.com](mailto:info@apexbt.com)

