


Animal Tissue Transparent Mounting Solution

Introduction

Animal Tissue Transparent Mounting Solution is an animal tissue clearing mounting medium that can be used in combination with the Animal Tissue Optical Clearing Kit.

Note

1. Storage conditions: Store at room temperature, valid for two years. It can also be stored at 4°C or -20°C.
2. This product is for scientific use only.



APExBIO Technology
www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com

