

Recombinant Human Melanoma Inhibitory Activity Protein

Information

Gene ID	8190
Accession #	Q16674
Alternate Names	
Source	Escherichia coli.
M.Wt	Approximately 12.1 kDa, a single non-glycosylated polypeptide chain containing 107 amino acids.
AA Sequence	GPMPKLADRK LCADQECSHP ISMAVALQDY MAPDCRFLTI HRGQVVYVFS KLKGRGRLFW GGSVQGDYYG DLAARLGYFP SSIVREDQTL KPGKVDVKTD KWDFYCQ
Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles - 12 months from date of receipt, -20 to -70 °C as supplied - 1 month, 2 to 8 °C under sterile conditions after reconstitution - 3 months, -20 to -70 °C under sterile conditions after reconstitution
Formulation	Lyophilized from a 0.2 μm filtered concentrated solution in PBS, pH 7.4, with 5 % Trehalose.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a concentration of 0.1-1.0 mg/ml. Stock solutions should be apportioned into working aliquots and stored at \leq -20 °C. Further dilutions should be made in appropriate buffered solutions.
Biological Activity	Fully biologically active when compared to standard. The ED as determined by a cell proliferation assay using human A375 cell line is less than 5 μ g/ml, corresponding to a specific activity of > 200 IU/mg.
Shipping Condition	Gel pack.
Handling	Centrifuge the vial prior to opening.

Components and Storage

Components	5µg	100µg	500µg
Recombinant Human Melanoma Inhibitory Activity Protein	5µg	100µg	500µg

Use a manual defrost freezer and avoid repeated freeze-thaw cycles

- 12 months from date of receipt, -20 to -70 °C as supplied
- 1 month, 2 to 8 °C under sterile conditions after reconstitution
- 3 months, -20 to -70 °C under sterile conditions after reconstitution

Quality Control

Purity	> 98 % by SDS-PAGE and HPLC analyses.	The state of the s
Endotoxin	Less than 0.1 EU/µg of rHuMIA as determine	d by LAL method.

Description

Reference





20





APExBIO Technology

www.apexbt.com

7505 Fannin street, Suite 410, Houston, TX 77054.

Tel: +1-832-696-8203 | Fax: +1-832-641-3177 | Email: info@apexbt.com